

INSTRUMED PVT LTD.,



Diagnostic Instruments

for ENT & Ophthalmic Examination

Diagnostic Instruments

for ENT & Ophthalmic Examination

Diagnostic Instruments Set

Full featured 2.5V set with 5 different instruments packed in hard plastic case

02-3050

Ophthalmoscope Head

- is equipped with 29 corrective lenses from -35 to 40 diopter

Otoscope Head

- Has fitting for insufflator
- Comes with removable 3x magnifying lens
- Ear speculum/cannula with internal diameter Ø 2, 3 & 4 mm, reusable, made from polycarbonate, autoclaveable at 121° C

Battery Handle

- Ensure connection to instrument heads with a bayonet locking system
- Require 2 "C" type batteries, not included
- Has knurled finish for better grip
- Has a rheostat to vary light intensity
- Spare lamp

Nasal Speculum

- Connects to the otoscope head and comes with manual screw for opening and closing of speculum

Bent Arm Illuminator

- Light the oral cavity and pharyngeal area
- Attaches to laryngeal mirrors No. 2 & 4 (included)
- Attaches to one metal tongue depressor (included)
- Attaches to one reusable adapter (included) that holds disposable tongue spatula



Otoscope

Packed in hard plastic case

02-3010

Otoscope Head

- Has fitting for insufflator
- Comes with removable 3x magnifying lens
- Ear speculum/cannula with internal diameter Ø 2, 3 & 4 mm, reusable, made from polycarbonate, autoclaveable at 121° C

Battery Handle

- Ensures connection to instrument heads with a bayonet locking system
- Requires 2 "C" type batteries, (not included)
- Has knurled finish for better grip
- Has a rheostat to vary light intensity
- Spare lamp





Otoscope / Ophthalmoscope Combination Set

Packed in hard plastic case

02-3030

Otoscope Head

- Has fitting for insufflator
- Comes with removable 3x magnifying lens
- Ear speculum/cannula with internal diameter Ø 2, 3 & 4 mm, reusable, made from polycarbonate, autoclaveable at 121° C

Ophthalmoscope Head

- is equipped with 29 corrective lenses from -35 to 40 diopter

Battery Handle

- Ensures connection to instrument heads with a bayonet locking system
- Requires 2 "C" type batteries, not included
- Has knurled finish for better grip
- Has a rheostat to vary light intensity
- Spare lamp



Ophthalmoscope

Packed in hard plastic case

02-3020

Ophthalmoscope Head

- is equipped with 29 corrective lenses from -35 to 40 diopter

Battery Handle

- Ensures connection to instrument heads with a bayonet locking system
- Requires 2 "C" type batteries, not included
- Has knurled finish for better grip
- Has a rheostat to vary light intensity
- Spare lamp



Diagnostic Instruments

for ENT & Ophthalmic Examination

Mini Otoscope / Ophthalmoscope Combination Set

Packed in hard plastic case

02-4030

- 2.5V LED illumination for true tissue color
- Wide angle swivel 2.5X viewing lens permits instrumentation
- 19 viewing lenses from -20 to +20 diopter for better resolution
- "AA" battery handle (batteries not included) with on/off power switch integrated into pocket clip
- 4 year speculum/cannula with internal diameter 2, 3, 4 & 5 mm



Mini Otoscope

02-4010

- 2.5V LED illumination for true tissue color
- Wide angle swivel 2.5X viewing lens permits instrumentation
- "AA" battery handle (batteries not included) with on/off power switch integrated into pocket clip
- 4 year speculum/cannula with internal diameter 2, 3, 4 & 5 mm
- Nylon pouch with belt loop



Mini Ophthalmoscope

02-4020

- 2.5V LED illumination for true tissue color
- 19 viewing lenses from -20 to +20 diopter for better resolution
- "AA" battery handle (batteries not included) with on/off power switch integrated into pocket clip
- Nylon pouch with belt loop





TAYLOR Percussion Hammer

Latex-free

	Weight	Length
02-6010	70g	20 cm/8"



BERLINER Percussion Hammer

Latex-free

	Weight	Length
02-6050	130g	20 cm/8"



DEJERINE Percussion Hammer

Latex-free

	Weight	Length
02-6020	113g	21 cm/8¼"



DEJERINE Percussion Hammer

Latex-free, with needle

	Weight	Length
02-6030	115g	21 cm/8¼"



TROEMNER Percussion Hammer

Latex-free

	Weight	Length
02-6040	210g	22 cm/8¾"




Diagnostic Instruments

Percussion Hammers

BUCK Percussion Hammer

Latex-free, with needle & brush

	Weight	Length
02-6060	110g	18 cm/7"



BABINSKY Percussion Hammer

Latex-free, head Ø 45 mm

	Weight	Length
02-6070	150g	21 cm/8½"



NATIONAL HOSPITAL Percussion Hammer

Latex-free, head Ø 50 mm

	Weight	Length
02-6080	95g	35 cm/14"





INSTRUMED PVT LTD,.

Diagnostic Instruments

for ENT & Ophthalmic Examination



Instrumed Pvt Ltd

65-B, S.I. Estate P.O. Box 615
Sialkot 51310, Pakistan

Tel +92 (52) 3551769

Fax +92 (52) 3559127

info@instrumed.com.pk

www.instrumed.com.pk

 **INSTRUMED PVT LTD,.**

www.instrumed.com.pk