

# PLASTIC SURGERY INSTRUMENTS



**Instrumed**  
Fine Surgical Instruments





Fiber Optic Retractors

1 - 6



Fiber Optic Light Guides

7



Retractors

8 - 15



Scissors

16 - 23



Dissectors

24 - 33



Markers

34 - 35



Calipers

36



Scalpels & Knives

37

## Plastic Surgery



Thumb & Tissue Forceps

38 - 41



Hemostatic Forceps

42



Speculas & Septum Forceps

43



Suction Tubes

44



Liposuction Cannulas

45 - 50



Index Alphabetical & Numerical

51 - 55



## Mammoplasty Retractors

### TEBBETS Style Retractors

with serrated tip

| With Light Guide | With Light Guide & Suction | Blades        |
|------------------|----------------------------|---------------|
| TR.120.1608      | TR.220.1608                | 16 mm x 8 cm  |
| TR.120.2409      | TR.220.2409                | 24 mm x 9 cm  |
| TR.120.2410      | TR.220.2410                | 24 mm x 10 cm |
| TR.120.2510      | TR.220.2510                | 25 mm x 10 cm |
| TR.120.2512      | TR.220.2512                | 25 mm x 12 cm |
| TR.120.2515      | TR.220.2515                | 25 mm x 15 cm |
| TR.120.3009      | TR.220.3009                | 30 mm x 9 cm  |
| TR.120.3012      | TR.220.3012                | 30 mm x 12 cm |
| TR.120.3015      | TR.220.3015                | 30 mm x 15 cm |
| TR.120.3016      | TR.220.3016                | 36 mm x 16 cm |
| TR.120.3019      | TR.220.3019                | 36 mm x 19 cm |

### FERRIERA Style Retractors

with smooth tip

| With Light Guide | With Light Guide & Suction | Blades        |
|------------------|----------------------------|---------------|
| FR.120.1407      | FR.220.1407                | 14 mm x 7 cm  |
| FR.120.1409      | FR.220.1409                | 14 mm x 9 cm  |
| FR.120.1608      | FR.220.1608                | 16 mm x 8 cm  |
| FR.120.1605      | FR.220.1605                | 16 mm x 15 cm |
| FR.120.2310      | FR.220.2310                | 23 mm x 10 cm |
| FR.120.2315      | FR.220.2315                | 23 mm x 15 cm |
| FR.120.2318      | FR.220.2318                | 23 mm x 18 cm |
| FR.120.2510      | FR.220.2510                | 25 mm x 10 cm |
| FR.120.2512      | FR.220.2512                | 25 mm x 12 cm |
| FR.120.2515      | FR.220.2515                | 25 mm x 15 cm |
| FR.120.2519      | FR.220.2519                | 25 mm x 19 cm |
| FR.120.3012      | FR.220.3012                | 30 mm x 12 cm |
| FR.120.3015      | FR.220.3015                | 30 mm x 15 cm |

### BIGGS Mammoplasty Retractors

| With Light Guide | With Light Guide & Suction | Retractor Only | Blades Width |
|------------------|----------------------------|----------------|--------------|
| FR.141.3018      | FR.241.3018                | FR.040.3018    | 30 mm Narrow |
| FR.141.5018      | FR.241.5018                | FR.040.5018    | 50 mm Wide   |



# Fiber Optic Retractors



**Instrumed**  
Fine Surgical Instruments

## Fiber Optic Retractors

Mutli Approach Retractor  
with smooth round tip

|             | Blade         |
|-------------|---------------|
| FR.020.2318 | 23 mm x 18 cm |

with Endoscopic Channel & Stopcock,  
to hold 10 mm telescopes.

|             | Blade         |
|-------------|---------------|
| FR.520.2518 | 25 mm x 18 cm |

Face-Lift Retractor

| Retractor Only | With Light Guide | Blade         |
|----------------|------------------|---------------|
| FR.045.2211    | FR.145.2211      | 22 mm x 11 cm |





## Fiber Optic Retractors

### MALINIAC Nasal Retractor

Length 4"/10 cm

| With Light Guide | Blade     |
|------------------|-----------|
| FR.146.1144      | 9 x 45 mm |



### CONVERSE Nasal Retractor

Length 3 1/2"/9 cm

| Retractor Only | With Light Guide | Blades     |
|----------------|------------------|------------|
| FR.060.1016    | FR.160.1016      | 10 x 20 mm |
| FR.060.1230    | FR.160.1230      | 12 x 30 mm |
| FR.060.1340    | FR.160.1340      | 12 x 40 mm |
| FR.060.1355    | FR.160.1355      | 12 x 50 mm |



### AUFRICHT Nasal Retractor

| Retractor Only | With Light Guide | Length       | Blades |
|----------------|------------------|--------------|--------|
| FR.050.0716    | FR.150.0716      | 6 1/4"/16 cm | 7 mm   |
| FR.050.1018    | FR.150.1018      | 7"/18 cm     | 10 mm  |



### WIEDER Tongue Depressors

Serrated Blade

| With Light Guide | With Light Guide & Suction | Depressore only | Blades      | Length        |
|------------------|----------------------------|-----------------|-------------|---------------|
| FR.165.3814      | FR.265.3814                | FR.065.3814     | 38 mm Small | 5 1/2" /14 cm |
| FR.165.4515      | FR.265.4515                | FR.065.4515     | 45 mm Large | 6" /15 cm     |



# Fiber Optic Retractors



**Instrumed**  
Fine Surgical Instruments

## Fiber Optic Retractors

### FREEMAN Flap Retractors

| Retractor Only | With Light Guide | Blades | Length   |
|----------------|------------------|--------|----------|
| FR.055.1320    | FR.155.1320      | 13 mm  | 7"/18 cm |
| FR.055.2520    | FR.155.2520      | 12 mm  | 8"/20 cm |
| FR.055.3520    | FR.155.3520      | 25 mm  | 8"/20 cm |

### HARRINGTON Retractors

12½"/32 cm

| Retractor Only | With Light Guide | With Light Guide & Suction | Blades      |
|----------------|------------------|----------------------------|-------------|
| FR.081.4013    | FR.181.4013      | FR.281.4013                | 40 x 127 mm |
| FR.081.6013    | FR.181.6013      | FR.281.6013                | 60 x 127 mm |

### DEAVER Retractors

| Retractor Only | With Light Guide | Blades    | Length    |
|----------------|------------------|-----------|-----------|
| FR.080.1918    | FR.180.1918      | ¾"/19 mm  | 7"/18 cm  |
|                | FR.180.1320      | ½"/13 mm  | 8"/20 cm  |
| FR.080.2220    | FR.180.2220      | ⅞"/22 mm  | 8"/20 cm  |
|                | FR.180.2520      | 1"/25 mm  | 8"/20 cm  |
| FR.080.2523    | FR.180.2523      | 1"/25 mm  | 9"/23 cm  |
| FR.080.2530    | FR.180.2530      | 1"/25 mm  | 12"/30 cm |
| FR.080.2533    |                  | 1"/25 mm  | 13"/33 cm |
| FR.080.2535    |                  | 1"/25 mm  | 14"/35 cm |
| FR.080.3830    | FR.180.3830      | 1½"/38 mm | 12"/30 cm |
| FR.080.5030    |                  | 2"/50 mm  | 12"/30 cm |
| FR.080.7530    | FR.180.7530      | 3"/75 mm  | 12"/30 cm |





## Extra / Intra Oral Retractors

### SIGMOID NOTCH Retractors

Lenght 8¾"/22 cm, Extra Oral Retractors

| Retractor Only | With Light Guide | Working Tip |
|----------------|------------------|-------------|
| FR.091.1625    | FR.191.1622      | 9.5/16 mm   |

### PETRI PTERYGOID

Lenght 9"/23 cm, Intra Oral Retractors

| Retractor Only | With Light Guide | Working Tip |
|----------------|------------------|-------------|
| FR.093.2527    | FR.193.2522      | 25.5 mm     |

### SHEA

Lenght 9¾"/25 cm, Multipurpose Retractor

| With Light Guide | Working Tip |
|------------------|-------------|
| FR.199.1425      | 14.5 mm     |



# Fiber Optic Retractors



**Instrumed**  
Fine Surgical Instruments

## Intra Oral Retractors

### BAUER

Length 8 3/4" / 21 cm, Intra Oral Retractors

| Retractor Only | With Light Guide | Working Tip                  |
|----------------|------------------|------------------------------|
| FR.094.1626    | FR.194.1626      | 26.5 cm / 16 mm, Left cheek  |
| FR.095.1626    | FR.195.1626      | 26.5 cm / 16 mm, Right cheek |

for Left Cheek

for Right Cheek

for Right Cheek

for Left Cheek

### LEAVASSEUR MERRILL Retractors

Length 9 3/4" / 25 cm, for Vertical Osteotomies

| Retractor Only | With Light Guide | Working Tip  |
|----------------|------------------|--------------|
| FR.097.2728    | FR.197.2725      | 12.5 / 27 mm |



### Light Guide Cables without Adaptors

| Black      | Gray       | Lenght          | Diameter |
|------------|------------|-----------------|----------|
| FC.124.BZZ | FC.124.GZZ | 6 Feet 183 cm   | Ø 3.5mm  |
| FC.144.BZZ | FC.144.GZZ | 7.5 Feet 229 cm | Ø 3.5mm  |
| FC.127.BZZ | FC.127.GZZ | 6 Feet 183 cm   | Ø 5.0mm  |
| FC.147.BZZ | FC.147.GZZ | 7.5 Feet 229 cm | Ø 5.0mm  |



### Light Source end Adaptors for Cables

| Adaptors   | Types                 |
|------------|-----------------------|
| FC.1A7.NON | Olympus               |
| FC.1A7.NSN | Storz                 |
| FC.1A7.NAN | ACMI British Standard |
| FC.1A7.NWN | Wolf                  |



### Instruments end Adaptors for Cables

| Adaptors   | Types                 |
|------------|-----------------------|
| FC.1A7.NNW | Wolf                  |
| FC.1A7.NNA | ACMI British Standard |
| FC.1A7.NWS | Storz                 |



# Retractors

## Breast Retractors



**Instrumed**  
Fine Surgical Instruments

### Breast Spatula with Handle

Length 12½"/31.5 cm, Semi malleable


|             | Working Tip     |
|-------------|-----------------|
| FR.056.2531 | 25 mm x 18.5 cm |
| FR.056.3831 | 38 mm x 18.5 cm |



### Breast Retractor

Length 12½"/31.5 cm, double-ended

|             | Blades                        |
|-------------|-------------------------------|
| FR.057.3831 | 22 mm x 55 mm , 38 mm x 70 mm |





## Rake Retractors

### FREEMAN Free-Lift Retractor

Lenght 7"/18 cm, with bent neck


|             | Blade | Prongs  |
|-------------|-------|---------|
| FR.072.3818 | 38 mm | 4 Sharp |

### FREEMAN Rake Retractors

Lenght 7"/18 cm

|             | Blade | Prongs               |
|-------------|-------|----------------------|
| FR.172.1918 | 19 mm | 4 blunt (2 recessed) |
| FR.173.1918 | 19 mm | 4 sharp (2 recessed) |
| FR.072.2218 | 22 mm | 4 sharp              |
| FR.172.2218 | 22 mm | 4 sharp (2 recessed) |
| FR.073.3818 | 38 mm | 4 sharp              |
| FR.172.3818 | 38 mm | 4 sharp (2 recessed) |

|             | Blade | Prongs               | Length        |
|-------------|-------|----------------------|---------------|
| FR.172.3810 | 38mm  | 4 sharp (2 recessed) | 4" (10 cm)    |
| FR.173.2510 | 25 mm | 4 sharp (2 recessed) | 4¼" (10.8 cm) |



### FREEMAN mini RAKE Retractor

Length 2½/5 cm, with thumb ring

|             | Blade | Prongs  |
|-------------|-------|---------|
| FR.072.1950 | 19 mm | 4 sharp |

### FREEMAN Rake Retractors

Lenght 4"/10 cm

|             | Blade |
|-------------|-------|
| FR.073.2210 | 22 mm |
| FR.073.3810 | 38 mm |

### GORNEY Rake Retractor

with spring neck

|             | Blade | Lenght      |
|-------------|-------|-------------|
| FR.074.3807 | 38 mm | 3"/7.5 cm   |
| FR.074.3808 | 38 mm | 4¼"/10.8 cm |



# Retractors

## Rake Retractors



**Instrumed**  
Fine Surgical Instruments

### MAXWELL Flap Retractors

Length 4"/10 cm

|             | Blades |
|-------------|--------|
| FR.075.2210 | 22 mm  |
| FR.075.3810 | 38 mm  |

### MILLARD Thimble Hooks

Blade 18 mm

| Single Hook | Double Hook | Prongs       |
|-------------|-------------|--------------|
| FR.076.1018 | FR.076.2018 | right handed |
| FR.176.1018 | FR.176.2018 | left handed  |

### Coronal Brow Lift Retractors

Length 7"/18 cm, 7 Prongs

|             | Blades | Grip hooks |
|-------------|--------|------------|
| FR.077.7018 | 70 mm  | front      |
| FR.177.7018 | 70 mm  | reverse    |

### Coronal Lift Retractor


Length 8 1/4"/21 cm, 6 Prongs

|             | Blade width |
|-------------|-------------|
| FR.078.1521 | 6"/15cm     |

### YANCOSKIE Abdominoplasty Retractor

Length 7 1/2"/19 cm, 6 prongs

|             | Blade width |
|-------------|-------------|
| FR.079.1519 | 6"/15cm     |





## Rake Retractors

### ANDERSON Bear Claw Retractors

|             |             | Length    | Prongs  |
|-------------|-------------|-----------|---------|
| FR.062.0307 | FR.162.0307 | 2¾"/7 cm  | 3 Sharp |
| FR.062.0511 | FR.162.0511 | 4½"/11 cm | 5 Sharp |

### CONVERSE Alar Retractor

Lenght 4¾"/12 cm, double ended

|             | Blade             |
|-------------|-------------------|
| FR.064.1044 | 10 mm/14mm, blunt |

### MARTINEZ KILNER Retractors

2 Sharp hooks, with thimble

|             | Length      | Tip            |
|-------------|-------------|----------------|
| FR.066.1008 | 3¼"/8 cm    | 10 mm (Swivel) |
| FR.066.7511 | 4½"/11.5 cm | 7.5 mm         |

### FREEMAN

double skin hook, sharp

|             | Length   |
|-------------|----------|
| FR.063.0210 | 4"/10 cm |

### RAGNELL-DAVIS Retractor

double ended

|             | Length    |
|-------------|-----------|
| RR.100.0514 | 5½"/14 cm |

### ROLLET Rake Retractor

Lenght 5"/13 cm

|             | Prongs  |
|-------------|---------|
| RR.120.0513 | 4 Sharp |

### MILLER Rake Retractors

Lenght 5¼"/13.5 cm

|             | Prongs  |
|-------------|---------|
| RR.130.0513 | 9 Sharp |
| RR.131.0513 | 9 Blunt |

### SENN Retractor

Lenght 6¼"/16 cm, Blade 7 x 8.5 mm

|             | Prongs  |
|-------------|---------|
| RR.150.0616 | 3 Sharp |



# Retractors



**Instrumed**  
Fine Surgical Instruments

## Lid Retractors / Skin Hooks / Nasal Hooks and Retractors

### KASDEN Retractor

double ended

|             | Length      |
|-------------|-------------|
| RR.155.0616 | 6½"/16.5 cm |



### DESMARRRES Lid Retractors

| Full Blade<br>5"/13 cm | Fenesrated<br>5"/13 cm | Full Blade<br>5½"/14 cm | Full Blade<br>6¼"/16 cm | Blade<br>size |
|------------------------|------------------------|-------------------------|-------------------------|---------------|
| RR.160.0813            |                        |                         | RR.166.0816             | 8 mm          |
| RR.160.1013            |                        | RR.164.1014             | RR.166.1016             | 10 mm         |
| RR.160.1213            | RR.162.1213            | RR.164.1214             | RR.166.1216             | 12 mm         |
| RR.160.1413            | RR.162.1413            | RR.164.1414             | RR.166.1416             | 14 mm         |
| RR.160.1613            | RR.162.1613            | RR.164.1614             | RR.166.1616             | 16 mm         |
| RR.160.1813            | RR.162.1813            |                         | RR.166.1816             | 18 mm         |

Fenestrated

Full



8 mm

10 mm

12 mm

14 mm


16 mm

18 mm

### JOSEPH Skin Hooks

Length 6¼"/16 cm

| Round Handle | Sharp / Prongs   |
|--------------|------------------|
| HK.010.0216  | 2 mm / 2 Prongs  |
| HK.010.0516  | 5 mm / 2 Prongs  |
| HK.010.0716  | 7 mm / 2 Prongs  |
| HK.010.1016  | 10 mm / 2 Prongs |



2 mm

5 mm

7 mm

10 mm

### COTTLE-NEIVERT Double Hook Retractor

double ended, with needle guide, blunt hooks

|             | Length    |
|-------------|-----------|
| RR.170.0821 | 8¼"/21 cm |



### FOMON Nasal Retractors

Length 6¼"/16 cm

|             |                        |
|-------------|------------------------|
| RR.181.0616 | single ball end prong  |
| RR.182.0616 | double ball end prongs |



### KLEINERT KUTZ Hooklets

Length 6¼"/16 cm

|             | Tips |
|-------------|------|
| HK.020.0316 | 3 mm |
| HK.020.0516 | 5 mm |



### GUTHRIE Hook Retractors

Length 5"/12.5 cm

|             | Fig | Tip    |
|-------------|-----|--------|
| HK.030.1012 |     | 1.0 mm |
| HK.030.1512 | 1   | 1.5 mm |
| HK.030.2012 | 2   | 2.0 mm |



### BARSKY Double Skin Hooks

|             | Length    | Width |
|-------------|-----------|-------|
| HK.350.0615 | 6"/15 cm  | 5 mm  |
| HK.390.0616 | 6¼"/16 cm | 9 mm  |



### TYRELL Skin Hook

Sharp

|             | Length   |
|-------------|----------|
| HK.410.0513 | 5"/13 cm |



### KILNER (TYRELL) Skin Hook

very delicate, with round handle

|             | Length   |
|-------------|----------|
| HK.420.0615 | 6"/15 cm |



### GILLIES Skin Hooks, delicate

Length 7"/18 cm

|             | Fig |
|-------------|-----|
| HK.431.0718 | 1   |
| HK.432.0718 | 2   |



# Retractors

## Facial / Nasal Retractors



**Instrumed**  
Fine Surgical Instruments

### BERNSTEIN Retractor

Length 5"/12.6 cm, double ended

|             | Blade     |
|-------------|-----------|
| FR.085.0846 | 8 x 46 mm |

### CLODIUS Orbital Retractor

double ended

|             |               |
|-------------|---------------|
| FR.082.1320 | 13 mm / 20 mm |
|-------------|---------------|

### WISE Orbital Retractor

double ended

|             | Length      |
|-------------|-------------|
| FR.083.0216 | 6½"/16.5 cm |

### CRONIN Cheek Retractor

double ended

|             | Length      |
|-------------|-------------|
| FR.087.0216 | 6½"/16.5 cm |

### FISCHL Shoe Horn Retractor

double ended

|             | Length    |
|-------------|-----------|
| FR.088.0216 | 6¼"/16 cm |





### SHEEN Retractor

Length 6"/15 cm

|             | Blade  |
|-------------|--------|
| FR.084.0115 | 6.5 mm |



### SHEEN Nasal Retractor

Length 6¼"/16 cm

|             | Blade |
|-------------|-------|
| FR.084.0116 | 10 mm |



### CRONIN Lip & Nostrill Retractor

double ended

|             | Length      |
|-------------|-------------|
| FR.087.0215 | 6¼"/15.5 cm |



### CAROLIN Finger Retractor

double ended

|             | Length    |
|-------------|-----------|
| FR.086.0211 | 4¼"/11 cm |



# Scissors



**Instrumed**  
Fine Surgical Instruments

## Rhytidectomy Scissors

### GORNEY Face-Lift Scissors

serrated blades

|          | TC          | SuperCut    | Length    |
|----------|-------------|-------------|-----------|
| Straight | SP.110.0719 | SP.120.0719 | 7½"/19 cm |
| Curved   | SP.111.0719 | SP.121.0719 | 7½"/19 cm |
| Straight | SP.110.0923 | SP.120.0923 | 9"/23 cm  |
| Curved   | SP.111.0719 | SP.121.0923 | 9"/23 cm  |

### FREEMAN Face-Lift Scissors

Length 7"/18 cm, serrated blades

|          | TC          | SuperCut    |
|----------|-------------|-------------|
| Straight | SP.130.0718 | SP.140.0718 |
| Curved   | SP.131.0718 | SP.141.0718 |

### FREEMAN-GORNEY Face-Lift Scissors

Length 7"/18 cm, serrated blades

|          | TC          | SuperCut    |
|----------|-------------|-------------|
| Straight | SP.150.0718 | SP.160.0718 |
| Curved   | SP.151.0718 |             |

### FREEMAN KAYE Face-Lift Scissors

Length 7"/18 cm, serrated blades

|          | TC          | SuperCut    |
|----------|-------------|-------------|
| Straight | SP.170.0718 | SP.180.0718 |
| Curved   | SP.171.0718 | SP.181.0718 |

### KAYE Face-Lift Scissors

serrated blades

|        | TC          | SuperCut    | Length      |
|--------|-------------|-------------|-------------|
| Curved | SP.211.0411 |             | 4½"/11.5 cm |
| Curved | SP.211.0615 | SP.221.0615 | 6"/15 cm    |
| Curved | SP.211.0718 | SP.221.0718 | 7"/18 cm    |
| Curved | SP.211.0923 | SP.221.0923 | 9"/23 cm    |





## Rhytidectomy and Nasal Scissors

### WILKINSON Face-Lift Scissors

Length 7½"/19 cm, Duckbill tip

|        | Standard    | TC          |
|--------|-------------|-------------|
| Curved | SP.230.0719 | SP.032.0719 |

### REES Face-Lift Scissors

one serrated blade

|          | TC          | SuperCut    | Length    |
|----------|-------------|-------------|-----------|
| Straight | SP.240.0617 | SP.250.0617 | 6¾"/17 cm |
| Curved   | SP.241.0617 | SP.251.0617 | 6¾"/17 cm |
| Straight | SP.240.0820 | SP.250.0820 | 8"/20 cm  |
| Curved   | SP.241.0820 | SP.251.0820 | 8"/20 cm  |

### ASTON Face-Lift Scissors

serrated blades

|          | TC          | Length    |
|----------|-------------|-----------|
| Straight | SP.260.0617 | 6¾"/17 cm |
| Curved   | SP.261.0617 | 6¾"/17 cm |
| Straight | SP.260.0820 | 8"/20 cm  |
| Curved   | SP.261.0820 | 8"/20 cm  |

### FOMON Lower Lateral Scissors

strongly curved

|  | SuperCut    | Length    |
|--|-------------|-----------|
|  | SC.010.0514 | 5½"/14 cm |

### FOMON Nasal Scissors

Length 5¾"/14.5 cm, one serrated blade

|        | TC          | SuperCut    |
|--------|-------------|-------------|
| Angled | SC.020.0514 | SC.022.0514 |



# Scissors



**Instrumed**  
Fine Surgical Instruments

## Nasal and Tenotomy Scissors

### FOMON Saber-back Dorsal scissors

flattened blunt tip, semi-sharp outer edges

|  | Angled 30°  | Angled 45°  | Length      |
|--|-------------|-------------|-------------|
|  | SC.035.0514 | SC.045.0514 | 5¾"/14.5 cm |



### KNIGHT Nasal Scissors

serrated

|        | TC          | SuperCut    | Length   |
|--------|-------------|-------------|----------|
| Angled | SC.040.0718 | SC.041.0718 | 7"/18 cm |

no image

### HEYMANN Nasal Scissors

serrated blades

|        | Standard    | Length   |
|--------|-------------|----------|
| Angled | SC.050.0718 | 7"/18 cm |



### COTTLE Nasal Scissors

one serrated blade

|        | TC          | SuperCut    | Length    |
|--------|-------------|-------------|-----------|
| Angled | SC.060.0616 | SC.062.0612 | 6¼"/16 cm |



### STEVENS Tenotomy Scissors

Length 4½"/11 cm

|          | TC          | SuperCut    |             |
|----------|-------------|-------------|-------------|
| Straight | SC.070.0411 | SC.080.0411 | Sharp/Sharp |
| Curved   | SC.071.0411 | SC.081.0411 | Sharp/Sharp |
| Straight | SC.072.0411 | SC.082.0411 | Blunt/Blunt |
| Curved   | SC.073.0411 | SC.083.0411 | Blunt/Blunt |





## Dissecting Scissors

### METZENBAUM Dissecting Scissors

|          | Standard    | TC          | SuperCut    | Length    |
|----------|-------------|-------------|-------------|-----------|
| Straight | SC.100.0514 | SC.110.0514 | SC.120.0514 | 5½"/14 cm |
| Curved   | SC.101.0514 | SC.111.0514 | SC.121.0514 | 5½"/14 cm |
| Straight | SC.100.0718 | SC.110.0718 | SC.120.0718 | 7"/18 cm  |
| Curved   | SC.101.0718 | SC.111.0718 | SC.121.0718 | 7"/18 cm  |
| Straight | SC.100.0923 | SC.110.0923 | SC.120.0923 | 9"/23 cm  |
| Curved   | SC.101.0923 | SC.111.0923 | SC.121.0923 | 9"/23 cm  |



### METZENBAUM Fino Scissors

Length 7"/18 cm

|          | TC          |
|----------|-------------|
| Straight | SC.090.0718 |
| Curved   | SC.091.0718 |



### LEXER-BABY Dissecting Scissors

Curved

| SuperCut    | Length   |
|-------------|----------|
| SC.095.0410 | 4"/10 cm |



### RAGNELL (KILNER) Dissecting Scissors

Curved

| TC          | SuperCut    | Length     |
|-------------|-------------|------------|
| SC.131.0512 | SC.133.0512 | 5"/12.7 cm |
|             | SC.133.0615 | 6"/15 cm   |



# Scissors

## Fine and Iris Scissors



**Instrumed**  
Fine Surgical Instruments

### IRIS Scissors

|          | Standard    | Length    |
|----------|-------------|-----------|
| Straight | SC.140.0309 | 3½"/9 cm  |
| Curved   | SC.141.0309 | 3½"/9 cm  |
| Straight | SC.140.0410 | 4⅞"/10 cm |
| Curved   | SC.141.0410 | 4⅞"/10 cm |

### IRIS Scissors

delicate, serrated, round screw joint

|          | SuperCut    | Length      |
|----------|-------------|-------------|
| Straight | SC.144.0411 | 4½"/11.5 cm |
| Curved   | SC.145.0411 | 4½"/11.5 cm |

### IRIS Scissors

Sharp/Sharp

|          | TC          | SuperCut    | Length      |
|----------|-------------|-------------|-------------|
| Straight | SC.148.0411 | SC.158.0411 | 4½"/11.5 cm |
| Curved   |             | SC.159.0411 | 4½"/11.5 cm |

### IRIS Ribbon Scissors

Sharp/Sharp

|          | Standard    | Length   |
|----------|-------------|----------|
| Straight | SC.150.0309 | 3½"/9 cm |
| Curved   | SC.151.0309 | 3½"/9 cm |
| Straight | SC.150.0410 | 4"/10 cm |
| Curved   | SC.151.0410 | 4"/10 cm |


### STRABISMUS Scissors

|          | TC          | SuperCut    | Length      |
|----------|-------------|-------------|-------------|
| Straight | SC.154.0411 | SC.164.0411 | 4½"/11.5 cm |
| Curved   | SC.155.0411 | SC.151.0410 | 4½"/11.5 cm |

### LOCKLIN Scissors

angled to side, sharp

|        | TC          | Length    |
|--------|-------------|-----------|
| Angled | SC.170.0616 | 6¼"/16 cm |





## Blepharoplasty Fine & Gum Scissors

### KAYE Blepharoplasty Scissors

ribbon handle, serrated

|        | Standard    | Length      |
|--------|-------------|-------------|
| Curved | SP.273.0411 | 4½"/11.5 cm |



### Blepharoplasty Scissors

flat tips, one serrated blade

|        | TC          | Length      |
|--------|-------------|-------------|
| Curved | SP.275.0411 | 4½"/11.5 cm |



### GRADLE Scissors

Sharp

|          | TC          | SuperCut    | Length     |
|----------|-------------|-------------|------------|
| Straight | SP.280.0309 | SP.290.0309 | 3¾"/9.5 cm |
| Curved   | SP.281.0309 | SP.291.0309 | 3¾"/9.5 cm |



### GRADLE Ribbon Scissors

Sharp

|          | Ribbon Handle | Length   |
|----------|---------------|----------|
| Straight | SP.282.0410   | 4"/10 cm |
| Curved   | SP.283.0410   | 4"/10 cm |



### WILMER Conjunctival and Utility Scissors

Angled on flat, Serrated cutting edge

|  | TC          |             | Length   |
|--|-------------|-------------|----------|
|  | SP.284.0410 | Blunt/Blunt | 4"/10 cm |
|  | SP.286.0411 | Sharp/Sharp | 4"/10 cm |



### GOLDMAN FOX Gum Scissors

|          | TC          | SuperCut    | Length   |
|----------|-------------|-------------|----------|
| Straight | SC.180.0513 | SC.190.0513 | 5"/13 cm |
| Curved   | SC.181.0513 | SC.191.0513 | 5"/13 cm |



# Scissors



**Instrumed**  
Fine Surgical Instruments

## Stitch & Wire Scissors, Micro Scissors

### O'BRIEN Stitch Scissors

fine angled tip, 15 mm cutting edge

|  | Standard    | Length    |
|--|-------------|-----------|
|  | SP.090.0410 | 4"/10 cm  |
|  | SP.290.0412 | 4¾"/12 cm |



### Universal Wire Scissors

one serrated blade

|          | TC          | Length    |
|----------|-------------|-----------|
| Serrated | SC.184.0412 | 4½"/12 cm |



### Micro Scissors, with round handle

|          |             | Length   |
|----------|-------------|----------|
| Straight | SM.010.0615 | 6"/15 cm |
| Curved   | SM.011.0615 | 6"/15 cm |
| Straight | SM.010.0718 | 7"/18 cm |
| Curved   | SM.011.0718 | 7"/18 cm |



### WESTCOTT Micro Scissors, with flat handle

Blade 10 mm, Sharp / Blunt

|        |             | Length    |
|--------|-------------|-----------|
| Curved | SM.021.0411 | 4¼"/11 cm |
| Curved | SM.031.0411 |           |





### GORNEY Scissors

Double Action

|        | TC          | Length     |
|--------|-------------|------------|
| Angled | SP.300.0820 | 8"/20.5 cm |



### CAPLAN Scissors

Double Action

|        | TC          | Length       |
|--------|-------------|--------------|
| Angled | SP.310.0821 | 8 1/4"/21 cm |



# Endo Dissectors

## Breast Dissectors



**Instrumed**  
Fine Surgical Instruments

### Spatulated Dissector

BD.045.0033      13"/33 cm



### AGRIS-DINGMAN Breast Dissectors

Length 14 1/4"/36 cm, graduated

BD.510.0036      Pair, (Set of 2)  
BD.010.0136      Left  
BD.010.0236      Right



### McCOLLUM-DINGMAN Submammary Dissector

BD.015.0031      12"/ 30 cm





## Breast Dissectors

### McCOLLUM-DINGMAN Breast Dissector

|             | Length                                 |
|-------------|----------------------------------------|
| BD.016.0041 | 16 <sup>1</sup> / <sub>4</sub> "/41 cm |



### REYNOLDS Transauxillary Breast Dissectors

Fenestrated

|             | Length    |
|-------------|-----------|
| BD.030.0033 | 13"/33 cm |
| BD.030.0040 | 16"/40 cm |



### Breast Dissectors

Angulated blades

|             | Length                                 |           |
|-------------|----------------------------------------|-----------|
| BD.050.0733 | 13"/33 cm                              | Rigid     |
| BD.050.0742 | 16 <sup>1</sup> / <sub>2</sub> "/42 cm | Rigid     |
| BD.350.0733 | 13"/33 cm                              | Malleable |
| BD.350.0742 | 16 <sup>1</sup> / <sub>2</sub> "/42 cm | Malleable |



# Endo Dissectors



**Instrumed**  
Fine Surgical Instruments

## Breast Dissectors

### Breast Dissectors

Spatulated blades

|             | Length     |           |
|-------------|------------|-----------|
| BD.051.3033 | 13"/33 cm  | Rigid     |
| BD.051.3042 | 16½"/42 cm | Rigid     |
| BD.351.3033 | 13"/33 cm  | Malleable |
| BD.351.3042 | 16½"/42 cm | Malleable |



### Breast Dissectors

Oval spatulated blades

|             | Length     |           |
|-------------|------------|-----------|
| BD.052.2533 | 13"/33 cm  | Rigid     |
| BD.052.2542 | 16½"/42 cm | Rigid     |
| BD.352.2533 | 13"/33 cm  | Malleable |
| BD.352.2542 | 16½"/42 cm | Malleable |



### EMORY Style Breast Dissector

Curved, blunt blade

|             | Length    |
|-------------|-----------|
| BD.055.1030 | 20"/30 cm |





## Forehead & Facelift Endo Dissectors, Scalp Elevators

### Scalp Elevator

Length 9½"/24 cm

|             | Shape           | Blade Width |
|-------------|-----------------|-------------|
| ED.111.1224 | Slightly curved | 12.5 mm     |



### Scalp Elevator

|             | Shape      | Blade Width |
|-------------|------------|-------------|
| ED.112.1224 | "S" Shaped | 12.5 mm     |



### Scalp Elevator

|             | Shape           | Blade Width |
|-------------|-----------------|-------------|
| ED.113.1524 | Slightly curved | 15 mm       |



### Scalp Elevator

|             | Shape           | Blade Width |
|-------------|-----------------|-------------|
| ED.114.1224 | Slightly curved | 12.5 mm     |



### Scalp Elevator

|             | Shape    | Blade Width |
|-------------|----------|-------------|
| ED.121.0724 | Straight | 7 mm        |



# Endo Dissectors




**Instrumed**  
Fine Surgical Instruments

## Forehead & Facelift Endo Dissectors, Scalp Elevators

### Scalp Elevator

Length 9½"/24 cm

|             | Shape          | Blade Width |
|-------------|----------------|-------------|
| ED.122.0724 | Quarter curved | 7 mm        |



### Scalp Elevator

|             | Shape        | Blade Width |
|-------------|--------------|-------------|
| ED.123.0724 | Fully curved | 7 mm        |



### Scalp Elevator

|             | Shape           | Blade Width |
|-------------|-----------------|-------------|
| ED.124.0728 | Slightly curved | 7 mm        |



### Scalp Elevator

Length 9¾"/24 cm, with Suction

|             | Shape           | Blade Width |
|-------------|-----------------|-------------|
| ED.622.0924 | Slightly curved | 7 mm        |



### Scalp Elevator

Length 10½"/26 cm

|             | Shape           | Blade Width |
|-------------|-----------------|-------------|
| ED.630.0826 | Slightly curved | 8.5 mm      |





## Frontotemporal, Frontoglabellar Dissectors

### Frontotemporal Dissector

Length 9"/22.5 cm

|             | Shape           | Blade Width |
|-------------|-----------------|-------------|
| ED.131.0522 | Slightly curved | 5 mm        |



### Frontotemporal Dissector

Length 9"/23 cm

|             | Shape           | Blade Width |
|-------------|-----------------|-------------|
| ED.132.0722 | Slightly curved | 7 mm        |



### Frontotemporal Dissector

Length 9"/23 cm

|             | Shape           | Blade Width |
|-------------|-----------------|-------------|
| ED.132.1022 | Slightly curved | 10 mm       |



### Frontoglabellar Dissector

Length 10½"/26.5 cm

|             | Shape      | Tip Size |
|-------------|------------|----------|
| ED.133.0726 | "S" Shaped | 7 mm     |



### Flap Dissector

Length 9¼"/23.5 cm

|             | Shape    | Tip Size |
|-------------|----------|----------|
| ED.134.0523 | Straight | 5 mm     |



# Endo Dissectors



**Instrumed**  
Fine Surgical Instruments

## Frontotemporal, Frontoglabellar Dissector, Periosteal Elevators

### Zygomatic Arch Elevator

Length 8¼"/20.5 cm

|             | Shape      | Tip Size |
|-------------|------------|----------|
| ED.135.0520 | Curved Tip | 5 mm     |



### Orbital Rim Dissector

Length 9"/22.5 cm

|             | Shape           | Blade Width |
|-------------|-----------------|-------------|
| ED.136.0522 | Slightly Curved | 7 mm        |



### Temporal T-Dissector

Length 9½"/23.5 cm

|             | Shape           | Tip Size |
|-------------|-----------------|----------|
| ED.140.1023 | Slightly Curved | 10 mm    |



### Periosteal Elevator

Length 10"/25 cm

|             | Shape     | Blade Width |
|-------------|-----------|-------------|
| ED.145.1225 | "S" Shape | 12 mm       |



### Periosteal Spreader

Size 9¾"/24 cm

|             | Shape     | Blade Width |
|-------------|-----------|-------------|
| ED.148.0824 | "S" Shape | 8 mm        |





## Periosteal Spreader

Size 9 $\frac{3}{4}$ "/24 cm


|             | Shape     | Blade Width |
|-------------|-----------|-------------|
| ED.148.1024 | "S" Shape | 10 mm       |



## Arcus Marginalis Dissector

Length 9 $\frac{3}{4}$ "/24 cm

|             | Shape      | Blade Width |
|-------------|------------|-------------|
| ED.149.1024 | "S" Shaped | 10 mm       |



## Midface Facial Dissector

Length 9 $\frac{3}{4}$ "/24 cm

|             | Shape           | Blade Width |
|-------------|-----------------|-------------|
| ED.150.1024 | Multiple curves | 12 mm       |



## Transoral Dissector

|             | Length                    |
|-------------|---------------------------|
| ED.151.0723 | 9 $\frac{3}{4}$ "/23.5 cm |



# Endo Dissectors

## Nerve Protector and Knife



**Instrumed**  
Fine Surgical Instruments

### Corrugated & Procerus Muscle Dissector

Length 7<sup>3</sup>/<sub>4</sub>"/19.5 cm, round handle

|             | Blade  |
|-------------|--------|
| ED.371.0519 | 4.7 mm |



### Nerve Protector

Length 6<sup>3</sup>/<sub>4</sub>"/17 cm

|            | Shape       | Blade  |
|------------|-------------|--------|
| ED373.0517 | Curved Left | 4.5 mm |



### Nerve Protector

Length 6<sup>3</sup>/<sub>4</sub>"/17 cm

|             | Shape        | Blade  |
|-------------|--------------|--------|
| ED.374.0517 | Curved Right | 4.5 mm |



### Nerve Cutting Knife

Length 6<sup>3</sup>/<sub>4</sub>"/17 cm

|             | Shape  | Blade |
|-------------|--------|-------|
| ED.375.0317 | Curved | 3 mm  |



### Knife Handle for, Plastic Surgery

Length 8<sup>1</sup>/<sub>4</sub>"/20.5 cm

ED.370.0020





### Nerve Hooks

Length 9<sup>3</sup>/<sub>4</sub>" / 24 cm

|             | Shaft     |
|-------------|-----------|
| ED.401.0124 | Rigid     |
| ED.401.0224 | Malleable |



### Nerve Hooks

Length 9<sup>1</sup>/<sub>2</sub>" / 24 cm

|             | Hooks |
|-------------|-------|
| ED.160.0124 | Left  |
| ED.160.0224 | Right |

left



right



### Nerve Dissector

Length 10<sup>1</sup>/<sub>2</sub>" / 26.5 cm

|             | Style   | Blade Width |
|-------------|---------|-------------|
| ED.155.0524 | "Cobra" | 4.7 mm      |



# Markers

## Breast Markers



**Instrumed**  
Fine Surgical Instruments

### GROSSMAN Areola Marker

|             | Size                          |
|-------------|-------------------------------|
| MA.410.0400 | 4 in 1 , 40 / 44 / 48 / 52 mm |

### FREEMAN Areola Markers

with inner markings, made from light weight aluminium

|             | Size - 2 in 1 |
|-------------|---------------|
| MA.220.3436 | 34 / 36 mm    |
| MA.220.3842 | 38 / 42 mm    |
| MA.220.4550 | 45 / 50 mm    |

### FREEMAN Areola Markers

stainless steel

|             | Size - 1 in 1 |
|-------------|---------------|
| MA.125.3800 | 38 mm         |
| MA.125.4200 | 42 mm         |
| MA.125.4500 | 45 mm         |
| MA.125.5000 | 50 mm         |

### FREEMAN Areola Markers

stainless steel, double sided

|             | Size - 2 in 1 |
|-------------|---------------|
| MA.225.3436 | 34 / 36 mm    |
| MA.225.3842 | 38 / 42 mm    |
| MA.225.4550 | 45 / 50 mm    |

### ACKERMANN Areola Markers


Length 4"/10 cm

|             | Size             |
|-------------|------------------|
| MA.030.0014 | 5½"/14 cm, small |
| MA.030.0016 | 6¼"/16 cm, large |

### McKISSOCK Areola Markers

Length 5½"/14 cm

|             | Size  |
|-------------|-------|
| MA.035.3800 | 38 mm |
| MA.035.4200 | 42 mm |
| MA.035.4500 | 45 mm |
| MA.035.5000 | 50 mm |





## Flap / Face lift Markers, Rulers and Grid

### PITANGUY Flap Grasping Forceps

|             | Length    |
|-------------|-----------|
| MA.955.0014 | 5½"/14 cm |

### D'ASSUMPCAO Face-Lift Markers

Length 6¼"/16 cm

| MA.960.0016 | without ratchet |
|-------------|-----------------|
| MA.960.0116 | with ratchet    |

### PITANGUY Facial Demarcator

|             | Length    |
|-------------|-----------|
| MA.956.0016 | 6¼"/16 cm |

### MARTEN Flap Markers

Length 5½"/14 cm

| MA.965.0014 | without ratchet |
|-------------|-----------------|
| MA.965.0114 | with ratchet    |

### Metal Rulers

with cm and inch graduations

|             | Measuring range |
|-------------|-----------------|
| MA.100.0815 | 6"/15 cm        |
| MA.100.0820 | 8"/20 cm        |

### SHEEN Grid

made from light weight aluminium

|             | Graduation |
|-------------|------------|
| MA.050.4070 | 70 x 40 mm |



# Calipers



**Instrumed**  
Fine Surgical Instruments

## Breast and Measuring Calipers

### McKISSOCK (COLLYER) Caliper

|             | Size      |
|-------------|-----------|
| CA.040.0030 | 12"/30 cm |



### HALPREN Breast Calipers

|             | Size  |
|-------------|-------|
| CA.010.3460 | 34 mm |
| CA.010.3860 | 38 mm |
| CA.010.4260 | 42 mm |



### JAMESON caliper

Scale, 0 - 80 mm

|             | Length          |
|-------------|-----------------|
| CA.100.0395 | 3 3/4" / 9.5 cm |



### THORPE (MARCHAC) caliper

Scale, 0 - 80 mm

|             | Length         |
|-------------|----------------|
| CA.200.0412 | 4 3/4" / 12 cm |



### CASTROVIEJO calipers

| Straight    | ScaleRang | Working lenght |
|-------------|-----------|----------------|
| CA.025.1590 | 0 - 15 mm | 3 1/2" / 9 cm  |
| CA.025.2009 | 0 - 20 mm | 3 1/2" / 9 cm  |





## Nasal Knives

### Scalpel handle

round

| Straight    | Handle | Size      |
|-------------|--------|-----------|
| SH.103.0514 | No. 3  | 5½"/14 cm |
| SH.104.0616 | No. 4  | 6¼"/16 cm |

### Scalpel handle

cm graduated, round

| Straight    | Handle | Size        |
|-------------|--------|-------------|
| SH.203.0514 | No. 3  | 5¾"/14.5 cm |

### JOSEPH Button End Knife

|             | Size     |
|-------------|----------|
| KN.100.0615 | 6"/15 cm |

### COTTLE Nasal Knife

|             | Size     |
|-------------|----------|
| KN.110.0615 | 6"/15 cm |

### FREER Mucosa and Knife

|             | Size     |
|-------------|----------|
| KN.120.0615 | 6"/15 cm |

### JOSEPH Nasal Saws

Length 7½"/19 cm

|             |               |
|-------------|---------------|
| KN.200.0719 | Straight      |
| KN.210.0719 | Left Bayonet  |
| KN.220.0719 | Right Bayonet |

### BALLENGER Swivel Knives

Length 7¾"/20 cm

| Straight    | Bayonet     | Blade Width |
|-------------|-------------|-------------|
| KN.303.0720 | KN.403.0720 | 3 mm        |
| KN.304.0720 | KN.404.0720 | 4 mm        |
| KN.305.0720 | KN.405.0720 | 5 mm        |
| KN.306.0720 | KN.306.0720 | 6 mm        |



# Forceps



**Instrumed**  
Fine Surgical Instruments

## Thumb & Tissue Forceps

### ADSON Dressing Forceps

Serrated tip

|             | Length                                  |
|-------------|-----------------------------------------|
| TF.010.0412 | 4 <sup>3</sup> / <sub>4</sub> " / 12 cm |
| TF.010.0615 | 6" / 15 cm                              |



### ADSON Dressing Forceps

Length 4<sup>3</sup>/<sub>4</sub>" / 12 cm, fenestrated lightweight handle

|             | Tip                     |
|-------------|-------------------------|
| TF.020.0412 | 2 mm serrated, straight |



### ADSON Tissue Forceps

| 1 x 2 Teeth | Tip  | Length                                  |
|-------------|------|-----------------------------------------|
| TF.030.0412 | 2 mm | 4 <sup>3</sup> / <sub>4</sub> " / 12 cm |
| TF.032.0412 |      | 4 <sup>3</sup> / <sub>4</sub> " / 12 cm |
| TF.030.0615 |      | 6" / 15 cm                              |



### ADSON Tissue Forceps

Length 4<sup>3</sup>/<sub>4</sub>" / 12 cm, fenestrated lightweight handle

|             | Shape                 |
|-------------|-----------------------|
| TF.040.0412 | 1 x 2 Teeth, straight |



### ADSON Dressing Forceps

TC

| Serrated    | Smooth      | Length                                  |
|-------------|-------------|-----------------------------------------|
| TF.050.0412 | TF.052.0412 | 4 <sup>3</sup> / <sub>4</sub> " / 12 cm |
| TF.050.0615 |             | 6" / 15 cm                              |



### ADSON Dressing Forceps

TC, Delicate

|  | Smooth      | Length                                  |
|--|-------------|-----------------------------------------|
|  | TF.056.0412 | 4 <sup>3</sup> / <sub>4</sub> " / 12 cm |





### ADSON Tissue Forceps

TC serrated tip, 1 x 2 Teeth

|             | Length                                  |
|-------------|-----------------------------------------|
| TF.060.0412 | 4 <sup>3</sup> / <sub>4</sub> " / 12 cm |
| TF.060.0615 | 6" / 15 cm                              |



### ADSON Micro Forceps

4<sup>3</sup>/<sub>4</sub>" / 12 cm, delicate

|             | Shape       |
|-------------|-------------|
| TF.070.0412 | Serrated    |
| TF.072.0412 | 1 x 2 Teeth |



### ADSON-BROWN Tissue Forceps

7 x 7 Teeth

|             | Length                                  |
|-------------|-----------------------------------------|
| TF.077.0412 | 4 <sup>3</sup> / <sub>4</sub> " / 12 cm |



### ADSON-BROWN Tissue Forceps

7 x 7 Teeth, fenestrated lightweight handle

|             | Length                                  |
|-------------|-----------------------------------------|
| TF.177.0412 | 4 <sup>3</sup> / <sub>4</sub> " / 12 cm |



### ADSON-BROWN Tissue Forceps

7 x 7 Teeth

| TC          | Length                                  |
|-------------|-----------------------------------------|
| TF.277.0412 | 4 <sup>3</sup> / <sub>4</sub> " / 12 cm |



### ADSON Tissue Forceps

Length 4<sup>3</sup>/<sub>4</sub>" / 12 cm, with tying platform

| Serrated    | Shape                |
|-------------|----------------------|
| TF.080.0412 |                      |
| TF.082.0412 | 1 x 2 Teeth          |
| TF.084.0412 | Delicate             |
| TF.086.0412 | 1 x 2 Teeth Delicate |



# Forceps



**Instrumed**  
Fine Surgical Instruments

## Blepharoplasty, Iris and Fine Forceps

### DANIEL-ADSON Forceps

Spoon / Cup shaped

|             | Length     |
|-------------|------------|
| TF.090.0512 | 5"/12.5 cm |

### ADSON Blepharoplasty Forceps

With delicate hooks

| TC          | Length   |
|-------------|----------|
| TP.010.0513 | 5"/13 cm |

### ADSON Blepharoplasty Forceps

With standard hooks and platform

| TC          | Length   |
|-------------|----------|
| TP.020.0513 | 5"/13 cm |

### CASTROVIEJO-BROWN Forceps

With multiple teeth

|             | Length      |
|-------------|-------------|
| TF.100.0411 | 4½"/11.5 cm |

### GREEN Fixation Forceps

Length 4"/10 cm

| Without lock | With lock   | Tip   |
|--------------|-------------|-------|
| TM.005.0410  | TM.007.0410 | 5 mm  |
| TM.010.0410  | TM.011.0410 | 10 mm |

### BISHOP-HARMON Eye Forceps

Length 3¼"/8.5 cm

|             | Tip                      |
|-------------|--------------------------|
| TM.028.0308 | 0.8 mm, 1 x 2 teeth      |
| TM.025.0308 | 0.5 mm, fine serrated    |
| TM.025.0312 | 0.5 mm, fine 1 x 2 teeth |





## Jeweller's Forceps, Polypus Forceps

### DESMARRES Chalazion Clamp

Length 3½"/9 cm, with fixation screw

|             | Size   | Tip        |
|-------------|--------|------------|
| TF.030.1117 | small  | 11 x 17 mm |
| TF.030.1223 | medium | 12 x 23 mm |
| TF.030.1728 | large  | 17 x 28 mm |



### JEWELLER'S Forceps

|              | Fig        | Length      | Tip                    |
|--------------|------------|-------------|------------------------|
| TM.0411.0001 | Fig. 1     | 4½"/11.5 cm | straight, standard     |
| TM.0411.0003 | Fig. 3     | 4½"/11.5 cm | straight, delicate     |
| TM.0411.0030 | Fig. 3C    | 4¼"/10.5 cm | straight, delicate     |
| TM.0411.0004 | Fig. 4     | 4¼"/11.0 cm | straight, very fine    |
| TM.0411.0005 | Fig. 5     | 4¼"/11.0 cm | straight, ultra fine   |
| TM.0411.0050 | Fig. 5A    | 4¼"/11.0 cm | angled, very fine      |
| TM.0411.0045 | Fig. 5/45° | 4¼"/11.0 cm | angled 45°, delicate   |
| TM.0411.0055 | Fig. 5/45° | 4¼"/11.0 cm | angled 45°             |
| TM.0411.0090 | Fig. 5/90° | 4¼"/11.0 cm | angled 90°, ultra fine |
| TM.0411.0007 | Fig. 7     | 4¼"/11.0 cm | curved, delicate       |



Fig. 1



Fig. 3



Fig. 3C



Fig. 4



Fig. 5



Fig. 7

### LUCAE Ear Polypus Forceps

Length 5½"/14 cm, bayonet

|             | Shape       |
|-------------|-------------|
| TF.050.0514 | Serrated    |
| TF.052.0514 | 1 x 2 Teeth |



### JANSEN-GRUENWALD Nasal Polypus Forceps

Bayonet

| Serrated    | Length      |
|-------------|-------------|
| TF.060.0616 | 6¼"/16 cm   |
| TF.060.0716 | 7½"/19 cm   |
| TF.060.0821 | 8½"/21.5 cm |



# Hemostatic Forceps



**Instrumed**  
Fine Surgical Instruments

## Tissue & Vulsellum Forceps, Towel Clamps

### ALLIS Tissue Forceps

Multiple teeth

| 4 x 5 Teeth | 5 x 6 Teeth |          |
|-------------|-------------|----------|
| HF.045.0513 |             | 5"/13 cm |
| HF.045.0615 | HF.056.0615 | 6"/15 cm |

### LAHEY Vulsellum Forceps

| 2 x 2 Prongs | 3 x 3 Prongs |          |
|--------------|--------------|----------|
| HF.022.0615  | HF.033.0615  | 6"/15 cm |
| HF.022.0820  | HF.033.0820  | 8"/20 cm |

### BACKHAUS Towel Clamps

|             | Length   |
|-------------|----------|
| HF.050.0309 | 3½"/9 cm |





## Speculas and Septum Forceps

**COTTLE Cartilage Crusher**  
With clip

NS.020.0000



**KILLIAN Nasal Speculas**  
Lenght 5"/13 cm

| Speculas Only | With Light Guide | Working Lenght |
|---------------|------------------|----------------|
| SP.010.3512   | SP.110.3512      | 35 mm          |
| SP.010.5012   | SP.110.5012      | 50 mm          |
| SP.010.7512   | SP.110.7512      | 75 mm          |
| SP.010.9012   | SP.110.9012      | 90 mm          |



**ASCH Septum Straightening Forceps**

NS.025.0023      9"/23 cm



**WALSHAM Septum Straightening Forceps**  
Lenght 9"/23 cm

|             |          |
|-------------|----------|
| NS.030.0023 | straight |
| NS.030.0123 | left     |
| NS.030.0223 | right    |



# Suction Tubes



**Instrumed**  
Fine Surgical Instruments

## Suction Tubes

### YANKAUER Suction Tube

complete with tube olive

| With Light Guide | Central bore | Working length |
|------------------|--------------|----------------|
| ST.120.1128      | Ø 2,0 mm     | 11¼"/28.5cm    |




### EICKEN Antrum Cannula

Length 6"/15 cm, Lure Lock

| Short curved | Long curved | Diameter |
|--------------|-------------|----------|
| ST.225.0615  | ST.325.0615 | Ø 2.5 mm |
| ST.230.0615  | ST.330.0615 | Ø 3.0 mm |
| ST.240.0615  | ST.340.0615 | Ø 4.0 mm |

short curved

long curved





### Liposuction Cannulas

Three holes, diagonal

| Cannulas    | Diameter | Working Length |
|-------------|----------|----------------|
| LS.101.0320 | Ø 3 mm   | 8"/20 cm       |
| LS.101.0330 | Ø 3 mm   | 12"/30 cm      |
| LS.101.0430 | Ø 4 mm   | 12"/30 cm      |
| LS.101.0530 | Ø 5 mm   | 12"/30 cm      |
| LS.101.0630 | Ø 6 mm   | 12"/30 cm      |



### Liposuction Cannulas

One central hole on upper side / two lateral holes below

| Cannulas    | Diameter | Working Length |
|-------------|----------|----------------|
| LS.102.0320 | Ø 3 mm   | 8"/20 cm       |
| LS.102.0330 | Ø 3 mm   | 12"/30 cm      |
| LS.102.0430 | Ø 4 mm   | 12"/30 cm      |
| LS.102.0530 | Ø 5 mm   | 12"/30 cm      |
| LS.102.0630 | Ø 6 mm   | 12"/30 cm      |



### Liposuction Cannulas

Two lateral holes on the upper side / one central hole below

| Cannulas    | Diameter | Working Length |
|-------------|----------|----------------|
| LS.103.0320 | Ø 3 mm   | 8"/20 cm       |
| LS.103.0330 | Ø 3 mm   | 12"/30 cm      |
| LS.103.0430 | Ø 4 mm   | 12"/30 cm      |
| LS.103.0530 | Ø 5 mm   | 12"/30 cm      |
| LS.103.0630 | Ø 6 mm   | 12"/30 cm      |



### Liposuction Cannulas

One long central hole on the upper side / two lateral holes below

| Cannulas    | Diameter | Working Length |
|-------------|----------|----------------|
| LS.105.0430 | Ø 4 mm   | 12"/30 cm      |
| LS.105.0530 | Ø 5 mm   | 12"/30 cm      |



### TOLEDO Liposuction Cannulas

Fork shaped, rounded tips, dissector between tips / one central hole

| Cannulas    | Diameter | Working Length |
|-------------|----------|----------------|
| LS.104.0205 | Ø 2 mm   | 2"/5 cm        |
| LS.104.0305 | Ø 3 mm   | 2"/5 cm        |
| LS.104.0405 | Ø 4 mm   | 2"/5 cm        |
| LS.104.0215 | Ø 2 mm   | 6"/15 cm       |
| LS.104.0315 | Ø 3 mm   | 6"/15 cm       |
| LS.104.0325 | Ø 3 mm   | 10"/25 cm      |
| LS.104.0405 | Ø 4 mm   | 10"/25 cm      |
| LS.104.0530 | Ø 5 mm   | 12"/30 cm      |



# Liposuction Cannulas



**Instrumed**  
Fine Surgical Instruments

## Infiltration Cannulas

### Liposuction Cannulas

Top tip with hole / one lateral hole

| Cannulas    | Diameter | Working Length |
|-------------|----------|----------------|
| LS.106.0320 | Ø 3 mm   | 8"/20 cm       |
| LS.106.0330 | Ø 3 mm   | 12"/30 cm      |



### Liposuction Cannulas

Half cut tip

| Cannulas    | Diameter | Working Length |
|-------------|----------|----------------|
| LS.107.0630 | Ø 6 mm   | 12"/30 cm      |
| LS.107.0635 | Ø 6 mm   | 14"/35 cm      |



### Liposuction Cannulas

10 holes, each hole Ø 1.0 mm, For infiltration

| Cannulas    | Diameter | Working Length |
|-------------|----------|----------------|
| LS.108.0225 | Ø 2 mm   | 10"/25 cm      |
| LS.108.0325 | Ø 3 mm   | 10"/25 cm      |



### MERCEDES Liposuction Cannulas

Three lateral holes

| Cannulas    | Diameter | Working Length |
|-------------|----------|----------------|
| LS.160.0215 | Ø 2 mm   | 6"/15 cm       |
| LS.160.0315 | Ø 3 mm   | 6"/15 cm       |
| LS.160.0325 | Ø 3 mm   | 10"/25 cm      |
| LS.160.0415 | Ø 4 mm   | 6"/15 cm       |
| LS.160.0225 | Ø 4 mm   | 10"/25 cm      |
| LS.160.0525 | Ø 5 mm   | 10"/25 cm      |
| LS.160.0530 | Ø 5 mm   | 12"/30 cm      |



### KOTZUR Liposuction Cannulas

9 holes

| Cannulas    | Diameter | Working Length |
|-------------|----------|----------------|
| LS.162.0330 | Ø 3 mm   | 12"/30 cm      |





### Liposuction Cannulas

One central hole, for fine facial contourings

| Handle Conn. | Luer-Lock Conn. | Diameter | Working Length |
|--------------|-----------------|----------|----------------|
| LS.164.1705  | LS.165.1705     | Ø 1.7 mm | 2"/5 cm        |
| LS.164.2005  | LS.165.2005     | Ø 2.0 mm | 2"/5 cm        |
| LS.164.0305  | LS.165.0305     | Ø 3.0 mm | 2"/5 cm        |
| LS.164.0210  | LS.165.0210     | Ø 2.0 mm | 4"/10 cm       |
| LS.164.1715  | LS.165.1715     | Ø 1.7 mm | 6"/15 cm       |
| LS.164.0215  | LS.165.0215     | Ø 2.0 mm | 6"/15 cm       |
| LS.164.0315  | LS.165.0315     | Ø 3.0 mm | 6"/15 cm       |



### Liposuction Cannulas

8 holes, for fine facial contourings

| Handle Conn. | Luer-Lock Conn. | Diameter | Dia.Holes | Working Length |
|--------------|-----------------|----------|-----------|----------------|
| LS.166.1710  |                 | Ø 1.7 mm | Ø 0.5 mm  | 4"/10 cm       |
| LS.166.0210  |                 | Ø 2.0 mm | Ø 1.0 mm  | 4"/10 cm       |
| LS.166.0215  | LS.167.0215     | Ø 2.0 mm | Ø 1.0 mm  | 6"/15 cm       |



### ENTNER Liposuction Cannulas

| Handle Conn. | Luer-Lock Conn. | Diameter | Dia.Holes | Holes | Working Length |
|--------------|-----------------|----------|-----------|-------|----------------|
| LS.170.0215  | LS.171.0215     | Ø 2.0 mm | Ø 1.0 mm  | 18    | 6"/15 cm       |
| LS.170.1215  | LS.171.1215     | Ø 2.5 mm | Ø 1.2 mm  | 18    | 6"/15 cm       |
| LS.170.1515  | LS.171.1515     | Ø 3.0 mm | Ø 1.5 mm  | 18    | 6"/15 cm       |
| LS.170.1220  | LS.171.1220     | Ø 2.5 mm | Ø 1.2 mm  | 22    | 8"/20 cm       |
| LS.170.1520  | LS.171.1520     | Ø 3.0 mm | Ø 1.5 mm  | 22    | 8"/20 cm       |
| LS.170.2020  | LS.171.2020     | Ø 4.0 mm | Ø 2.0 mm  | 22    | 8"/20 cm       |
| LS.170.1215  | LS.171.1215     | Ø 2.5 mm | Ø 1.2 mm  | 26    | 10"/15 cm      |
| LS.170.1530  | LS.171.1530     | Ø 3.0 mm | Ø 1.5 mm  | 30    | 12"/30 cm      |
| LS.170.0230  | LS.171.0230     | Ø 4.0 mm | Ø 2.0 mm  | 30    | 12"/30 cm      |



# Liposuction Cannulas




**Instrumed**  
Fine Surgical Instruments

## Liposuction Cannulas

### ENTNER Liposuction Cannulas

22 holes, each hole Ø 1.5 mm


| Handle Conn. | Luer-Lock Conn. | Diameter | Working Length |
|--------------|-----------------|----------|----------------|
| LS.172.0320  | LS.173.0320     | Ø 3.0 mm | 8"/20 cm       |



### ENTNER Liposuction Cannulas

Bayonet shaped

| Handle Conn. | Luer-Lock Conn. | Diameter | Dia.Holes | Holes | Working Length |
|--------------|-----------------|----------|-----------|-------|----------------|
| LS.174.0320  | LS.175.0320     | Ø 3.0 mm | Ø 1.5 mm  | 22    | 8"/20 cm       |
| LS.174.0420  | LS.175.0420     | Ø 4.0 mm | Ø 2.0 mm  | 22    | 8"/20 cm       |
| LS.174.0330  | LS.175.0330     | Ø 3.0 mm | Ø 1.5 mm  | 30    | 12"/30 cm      |
| LS.174.0430  | LS.175.0430     | Ø 4.0 mm | Ø 2.0 mm  | 30    | 12"/30 cm      |





### Handles for Cannulas

With connection thread for cannulas hose connection

| Cannulas    | Diameter | Length    |
|-------------|----------|-----------|
| LS.120.1300 | Ø 13 mm  | 5¼"/13 cm |




| Cannulas    | Diameter | Length    |
|-------------|----------|-----------|
| LS.130.1300 | Ø 13 mm  | 5¼"/13 cm |




### Handles for Cannulas

With connection thread for cannulas hose connection

| Cannulas    | Diameter | Length    |
|-------------|----------|-----------|
| LS.121.1300 | Ø 13 mm  | 5½"/14 cm |



| Cannulas    | Diameter | Length    |
|-------------|----------|-----------|
| LS.131.1300 | Ø 13 mm  | 5½"/14 cm |



# Liposuction Cannulas



**Instrumed**  
Fine Surgical Instruments

## Thumb Plate Cannulas

### GASPAROTTI Liposuction Cannulas

One central hole, with attached Handle

| Straight    | Diameter | Working Length |
|-------------|----------|----------------|
| LS.151.0315 | Ø 3 mm   | 6"/15 cm       |
| LS.151.0615 | Ø 6 mm   | 6"/15 cm       |
| LS.151.0425 | Ø 4 mm   | 10"/25 cm      |
| LS.151.0625 | Ø 6 mm   | 10"/25 cm      |
| LS.151.0930 | Ø 9 mm   | 12"/30 cm      |



### GASPAROTTI Liposuction Cannulas

Two central hole, with attached Handle

| Straight    | Diameter | Working Length |
|-------------|----------|----------------|
| LS.153.0415 | Ø 4 mm   | 6"/15 cm       |
| LS.153.0625 | Ø 6 mm   | 10"/25 cm      |
| LS.153.0630 | Ø 6 mm   | 12"/30 cm      |



### GASPAROTTI Liposuction Cannulas

One central hole, with attached Handle

| Angled      | Diameter | Working Length |
|-------------|----------|----------------|
| LS.152.0315 | Ø 3 mm   | 6"/15 cm       |
| LS.152.0615 | Ø 6 mm   | 6"/15 cm       |
| LS.152.0625 | Ø 6 mm   | 10"/25 cm      |





## Index Alphabetical

### A

|                                   |       |
|-----------------------------------|-------|
| ACKERMANN Areola Markers          | 34    |
| ADSON Blepharoplasty Forceps      | 40    |
| ADSON-BROWN Tissue Forceps        | 39    |
| ADSON Dressing Forceps            | 38    |
| ADSON Micro Forceps               | 39    |
| ASCH Septum Straightening Forceps | 43    |
| ADSON Tissue Forceps              | 38-39 |
| ADAPTORS for Cables               | 7     |
| AGRIS-DINGMAN Breast Dissectors   | 24    |
| ALLIS Tissue Forceps              | 42    |
| ANDERSON Bear Claw Retractors     | 11    |
| ARCUS Marginalis Dissectors       | 31    |
| ASTON Face-Lift Scissors          | 17    |
| AUFRICHT Nasal Retractor          | 3     |

### B

|                                 |       |
|---------------------------------|-------|
| BACKHAUS Towel Clamps           | 42    |
| BALLENGER Swivel Knives         | 37    |
| BARSKY Double Skin Hooks        | 13    |
| BAUER Retractor                 | 6     |
| BERNSTEIN Retractors            | 14    |
| BISHOP-HARMON Eye Forceps       | 40    |
| BIGGS Mammoplasty Retractors    | 1     |
| Blepharoplasty Scissors         | 21    |
| Breast Spatula                  | 8     |
| Breast Retractor (Double-ended) | 8     |
| Breast Dissectors               | 24-26 |

### C

|                                        |    |
|----------------------------------------|----|
| Calipers                               | 36 |
| CAPLAN Scissors                        | 23 |
| CAROLIN Finger Retractor               | 15 |
| CASTROVIEJO-BROWN Forceps              | 40 |
| CASTROVIEJO Caliper                    | 36 |
| CLODIUS Orbital Retractor              | 14 |
| CONVERSE Nasal Retractor               | 3  |
| CORONAL Brow Lift Retractor            | 10 |
| CORONAL Lift Retractor                 | 10 |
| CONVERSE Alar Retractor                | 11 |
| COTTLE Cartilage Crusher               | 43 |
| COTTLE Nasal Knife                     | 37 |
| COTTLE-NEIVERT Double Hook Retractor   | 12 |
| COTTLE Nasal Scissors                  | 18 |
| Corrugated & Procerus Muscle Dissector | 32 |
| CRONIN Lips & Nostril Retractor        | 15 |
| CRONIN Cheek Retractor                 | 14 |

### D

|                              |    |
|------------------------------|----|
| DANIEL-ADSON Forceps         | 40 |
| D'ASSUMPCAO Face-Lift Marker | 35 |
| DEAVER Retractors            | 4  |
| DESMARRES Chalazion Clamp    | 41 |
| DESMARRES Lid Retractors     | 12 |

### E

|                             |       |
|-----------------------------|-------|
| ENDO DISSECTORS             | 24-33 |
| Extra Oral Retractors       | 6     |
| EICKEN Antrum Cannula       | 44    |
| ENTNER Liposuction Cannulas | 47-48 |

### F

|                           |   |
|---------------------------|---|
| FACE-LIFT Retractor       | 2 |
| FERRIERA Style Retractors | 1 |
| Fiberoptic Light Guide    | 7 |

|                                   |     |
|-----------------------------------|-----|
| Fiber Optic Retractors            | 1-6 |
| FISCHL Shoe Horn Retractor        | 14  |
| Flap Dissector                    | 29  |
| FOMON Lower Lateral Scissors      | 17  |
| FOMON Nasal Retractors            | 12  |
| FOMON Nasal Scissors              | 17  |
| FOMON Saber-back Dorsal scissors  | 18  |
| FREEMAN Areola Markers            | 34  |
| FREEMAN Flap Retractors           | 4   |
| FREEMAN Free-Lift Retractor       | 9   |
| FREEMAN Face-Lift Scissors        | 16  |
| FREEMAN-GORNEY Face-Lift Scissors | 16  |
| FREEMAN Hook Retractor            | 14  |
| FREEMAN KAYE Face-Lift Scissors   | 16  |
| FREEMAN Rake Retractor            | 9   |
| FREER Mucosa and Knife            | 37  |
| Frontotemporal Dissectors         | 29  |
| Frontoglabellar Dissectors        | 29  |

### G

|                                 |    |
|---------------------------------|----|
| GASPAROTTI Liposuction Cannulas | 50 |
| GILLIES Skin Hooks              | 13 |
| GOLDMAN FOX Gum Scissors        | 21 |
| GORNEY Face-Lift Scissors       | 16 |
| GORNEY Rake Retractor           | 9  |
| GORNEY Scissors                 | 23 |
| GRADLE Ribbon Scissors          | 21 |
| GRADLE Scissors                 | 21 |
| GREEN Fixation Forceps          | 40 |
| GROSSMAN Areola Marker          | 34 |
| GUTHRIE Hook Retractors         | 13 |

### H

|                        |    |
|------------------------|----|
| Handles for Cannulas   | 49 |
| HALPREN Breast Caliper | 36 |
| HARRINGTON Retractors  | 4  |
| HEYMANN Nasal Scissors | 18 |

### I

|                          |     |
|--------------------------|-----|
| Instruments End Adaptors | 7   |
| Intra Oral Retractors    | 6-5 |
| IRIS Ribbon Scissors     | 20  |
| IRIS Scissors            | 20  |

### J

|                                      |    |
|--------------------------------------|----|
| JAMESON caliper                      | 36 |
| JANSEN-GRUENWALD Nasal Polypus Force | 41 |
| JEWELLER'S Forceps                   | 41 |
| JOSEPH Button End Knife              | 37 |
| JOSEPH Nasal Saws                    | 37 |
| JOSEPH Skin Hooks                    | 12 |

### K

|                                  |    |
|----------------------------------|----|
| KASDEN Retractor                 | 12 |
| KAYE Blepharoplasty Scissors     | 21 |
| KAYE Face-Lift Scissors          | 16 |
| KILLIAN Nasal Speculas           | 43 |
| KILNER (TYRELL) Skin Hook        | 13 |
| KLEINERT KUTZ Hooklets           | 13 |
| KNIFE Handle for Plastic Surgery | 32 |
| KNIGHT Nasal Scissors            | 18 |
| KOTZUR Liposuction Cannulas      | 46 |

### L

|                         |    |
|-------------------------|----|
| LAHEY Vulsellum Forceps | 42 |
|-------------------------|----|

|                                |       |
|--------------------------------|-------|
| Light Source End Adaptors      | 7     |
| LEVASSEUR Merrill Retractor    | 6     |
| LEXER-BABY Dissecting Scissors | 19    |
| Light Guide Cables             | 7     |
| Liposuction Cannulas           | 45-48 |
| LOCKLIN Scissors               | 20    |
| LUCAE Ear Polypus Forceps      | 41    |

### M

|                                       |       |
|---------------------------------------|-------|
| MALINIAC Nasal Retractor              | 3     |
| Markers                               | 34,35 |
| MARTINEZ Kilner Retractor             | 11    |
| MARTEN Flap Marker                    | 35    |
| MAXWELL Flap Retractor                | 10    |
| MCCOLLUM-Dingman Submammary Dissector | 24    |
| MCCOLLUM-Dingman Breast Dissector     | 25    |
| MCKISSOCK Areola Markers              | 34    |
| MCKISSOCK (Collyer) Caliper           | 36    |
| Metal Rulers                          | 35    |
| MERCEDES Liposuction Cannulas         | 46    |
| METZENBAUM Dissecting Scissors        | 19    |
| METZENBAUM Fino Scissors              | 19    |
| Micro Scissors                        | 22    |
| Midface Facial Dissectors             | 31    |
| MILLER Rake Retractors                | 11    |
| Millard Thimble Hooks                 | 10    |
| Multi Approach Retractor              | 2     |

### N

|                     |    |
|---------------------|----|
| Nerve Protector     | 32 |
| Nerve Cutting Knife | 32 |
| Nerve Hook          | 33 |
| Nerve Dissector     | 33 |

### O

|                         |    |
|-------------------------|----|
| O'BRIEN Stitch Scissors | 22 |
| ORBITAL Rim Dissectors  | 30 |

### P

|                                |    |
|--------------------------------|----|
| PETRI PTERYGOID Retractor      | 5  |
| Periosteal Elevators           | 30 |
| Periosteal Spreaders           | 31 |
| PITANGUY Flap Grasping Forceps | 35 |
| PITANGUY Facial Demarcator     | 35 |

### R

|                                         |      |
|-----------------------------------------|------|
| RAGNELL-DAVIS Retractor                 | 11   |
| RAGNELL (KILNER) Dissecting Scissors    | 19   |
| Retractors                              | 1-15 |
| REES Face-Lift Scissors                 | 17   |
| REYNOLDS Transaxillary Breast Dissector | 25   |
| ROLLET Rake Retractor                   | 11   |

### S

|                           |       |
|---------------------------|-------|
| Scalp Elevators           | 27-28 |
| Scalpel handle            | 37    |
| Scissors                  | 16-23 |
| SENN Retractor            | 11    |
| SHEA Retractors           | 5     |
| SHEEN Grid                | 35    |
| SHEEN Nasal Retractor     | 15    |
| SHEEN Retractors          | 15    |
| Spatulated Dissector      | 24    |
| SIGMOID NOTCH Retractors  | 5     |
| STEVENS Tenotomy Scissors | 18    |



## Index Alphabetical

STRABISMUS Scissors .....20

### T

TEBBETS Style Retractors .....1

Temporal T-Dissector .....30

THORPE (MARCHAC) caliper .....36

TOLEDO Liposuction Cannulas .....46

Transoral Dissectors .....31

TYRELL Skin Hook .....13

### U

Universal Wire Scissors .....22

### W

WALSHAM Septum Straightening Forceps 43

WESTCOTT Micro Scissors .....22

WIEDER Tongue Depressors .....3

WILMER Conjunctival and Utility Scissors 21

WILKINSON Face-Lift Scissors .....17

WISE Orbital Retractor .....14

### Y

YANCOSKIE Abdominoplasty Retractor ...10

YANKAUER Suction Tube .....44

### Z

ZYGOMATIC Arch Elevator .....30



## Index Numerical

| Item               | Page      | Item        | Page | Item        | Page | Item        | Page | Item        | Page |
|--------------------|-----------|-------------|------|-------------|------|-------------|------|-------------|------|
| <b>BD</b>          |           | ED.375.0317 | 32   | FR.080.2220 | 4    | FR.180.1320 | 4    | KN.100.0615 | 37   |
| <b>BD.010.0136</b> | <b>24</b> | ED.401.0124 | 33   | FR.080.2523 | 4    | FR.180.2220 | 4    | KN.110.0615 | 37   |
| BD.010.0236        | 24        | ED.401.0224 | 33   | FR.080.2530 | 4    | FR.180.2520 | 4    | KN.120.0615 | 37   |
| BD.015.0031        | 24        | ED.622.0924 | 28   | FR.080.2533 | 4    | FR.180.2523 | 4    | KN.200.0719 | 37   |
| BD.016.0041        | 25        | ED.630.0826 | 28   | FR.080.2535 | 4    | FR.180.2530 | 4    | KN.210.0719 | 37   |
| BD.030.0033        | 25        | <b>FC</b>   |      | FR.080.3830 | 4    | FR.180.3830 | 4    | KN.220.0719 | 37   |
| BD.030.0040        | 25        | FC.124.BZZ  | 7    | FR.080.5030 | 4    | FR.180.7530 | 4    | KN.303.0720 | 37   |
| BD.045.0033        | 24        | FC.144.BZZ  | 7    | FR.080.7530 | 4    | FR.181.4013 | 4    | KN.304.0720 | 37   |
| BD.050.0733        | 25        | FC.127.BZZ  | 7    | FR.081.4013 | 4    | FR.181.6013 | 4    | KN.305.0720 | 37   |
| BD.050.0742        | 25        | FC.147.BZZ  | 7    | FR.081.6013 | 4    | FR.191.1622 | 5    | KN.306.0720 | 37   |
| BD.051.3033        | 26        | FC.124.GZZ  | 7    | FR.082.1320 | 14   | FR.193.2522 | 5    | KN.403.0720 | 37   |
| BD.051.3042        | 26        | FC.144.GZZ  | 7    | FR.083.0216 | 14   | FR.194.1626 | 6    | KN.404.0720 | 37   |
| BD.052.2533        | 26        | FC.127.GZZ  | 7    | FR.084.0115 | 15   | FR.195.1626 | 6    | KN.405.0720 | 37   |
| BD.052.2542        | 26        | FC.147.GZZ  | 7    | FR.084.0116 | 15   | FR.197.2725 | 6    | KN.406.0720 | 37   |
| BD.055.1030        | 26        | FC.1A7.NON  | 7    | FR.085.0846 | 14   | FR.199.1425 | 5    | <b>LS</b>   |      |
| BD.350.0733        | 25        | FC.1A7.NSN  | 7    | FR.086.0211 | 15   | FR.220.1407 | 1    | LS.101.0320 | 45   |
| BD.350.0742        | 25        | FC.1A7.NAN  | 7    | FR.087.0215 | 15   | FR.220.1409 | 1    | LS.101.0330 | 45   |
| BD.351.3033        | 26        | FC.1A7.NWN  | 7    | FR.087.0216 | 14   | FR.220.1605 | 1    | LS.101.0430 | 45   |
| BD.351.3042        | 26        | FC.1A7.NNW  | 7    | FR.088.0216 | 14   | FR.220.1608 | 1    | LS.101.0530 | 45   |
| BD.352.2533        | 26        | FC.1A7.NNA  | 7    | FR.091.1625 | 5    | FR.220.2310 | 1    | LS.101.0630 | 45   |
| BD.352.2542        | 26        | FC.1A7.NWS  | 7    | FR.093.2527 | 5    | FR.220.2315 | 1    | LS.102.0320 | 45   |
| BD.510.0036        | 24        | <b>FR</b>   |      | FR.094.1626 | 6    | FR.220.2318 | 1    | LS.102.0330 | 45   |
| <b>CA</b>          |           | FR.020.2318 | 2    | FR.095.1626 | 6    | FR.220.2510 | 1    | LS.102.0430 | 45   |
| CA.010.3460        | 36        | FR.040.3018 | 1    | FR.097.2728 | 6    | FR.220.2512 | 1    | LS.102.0430 | 45   |
| CA.010.3860        | 36        | FR.040.5018 | 1    | FR.120.1407 | 1    | FR.220.2515 | 1    | LS.102.0530 | 45   |
| CA.010.4260        | 36        | FR.045.2211 | 2    | FR.120.1409 | 1    | FR.220.2519 | 1    | LS.102.0630 | 45   |
| CA.025.1590        | 36        | FR.050.0716 | 3    | FR.120.1605 | 1    | FR.220.3012 | 1    | LS.103.0320 | 45   |
| CA.025.2009        | 36        | FR.050.1018 | 3    | FR.120.1608 | 1    | FR.220.3015 | 1    | LS.103.0330 | 45   |
| CA.040.0030        | 36        | FR.055.1320 | 4    | FR.120.2310 | 1    | FR.241.3018 | 1    | LS.103.0430 | 45   |
| CA.100.0395        | 36        | FR.055.2520 | 4    | FR.120.2315 | 1    | FR.241.5018 | 1    | LS.103.0530 | 45   |
| CA.200.0412        | 36        | FR.055.3520 | 4    | FR.120.2318 | 1    | FR.265.3814 | 3    | LS.103.0630 | 45   |
| <b>ED</b>          |           | FR.056.2531 | 8    | FR.120.2510 | 1    | FR.265.4515 | 3    | LS.104.0205 | 45   |
| ED.111.1224        | 27        | FR.056.3831 | 8    | FR.120.2512 | 1    | FR.281.4013 | 4    | LS.104.0305 | 45   |
| ED.112.1224        | 27        | FR.057.3831 | 8    | FR.120.2515 | 1    | FR.281.6013 | 4    | LS.104.0405 | 45   |
| ED.113.1524        | 27        | FR.060.1016 | 3    | FR.120.2519 | 1    | FR.520.2518 | 2    | LS.104.0215 | 45   |
| ED.114.1224        | 27        | FR.060.1230 | 3    | FR.120.3012 | 1    | <b>HF</b>   |      | LS.104.0315 | 45   |
| ED.121.0724        | 27        | FR.060.1340 | 3    | FR.120.3015 | 1    | HF.022.0615 | 42   | LS.104.0325 | 45   |
| ED.122.0724        | 28        | FR.060.1355 | 3    | FR.141.3018 | 1    | HF.022.0820 | 42   | LS.104.0405 | 45   |
| ED.123.0724        | 28        | FR.062.0307 | 11   | FR.141.5018 | 1    | HF.033.0615 | 42   | LS.104.0530 | 45   |
| ED.124.0728        | 28        | FR.062.0511 | 11   | FR.145.2211 | 2    | HF.033.0820 | 42   | LS.105.0430 | 45   |
| ED.131.0522        | 29        | FR.063.0210 | 11   | FR.146.1144 | 3    | HF.045.0513 | 42   | LS.105.0530 | 45   |
| ED.132.0722        | 29        | FR.064.1044 | 11   | FR.150.0716 | 3    | HF.045.0615 | 42   | LS.106.0320 | 46   |
| ED.132.1022        | 29        | FR.065.3814 | 3    | FR.150.1018 | 3    | HF.050.0309 | 42   | LS.106.0330 | 46   |
| ED.133.0726        | 29        | FR.065.4515 | 3    | FR.155.1320 | 4    | HF.056.0615 | 42   | LS.107.0630 | 46   |
| ED.134.0523        | 29        | FR.066.1008 | 11   | FR.155.2520 | 4    | <b>HK</b>   |      | LS.107.0635 | 46   |
| ED.135.0520        | 30        | FR.066.7511 | 11   | FR.160.1016 | 3    | HK.010.0216 | 12   | LS.108.0225 | 46   |
| ED.136.0522        | 30        | FR.072.2218 | 9    | FR.160.1230 | 3    | HK.010.0516 | 12   | LS.108.0325 | 46   |
| ED.140.1023        | 30        | FR.072.3818 | 9    | FR.160.1340 | 3    | HK.010.0716 | 12   | LS.120.1300 | 49   |
| ED.145.1225        | 30        | FR.072.1950 | 9    | FR.160.1355 | 3    | HK.010.1016 | 12   | LS.121.1300 | 49   |
| ED.148.0824        | 30        | FR.073.2210 | 9    | FR.162.0307 | 11   | HK.020.0316 | 13   | LS.130.1300 | 49   |
| ED.148.1024        | 31        | FR.073.3810 | 9    | FR.162.0511 | 11   | HK.020.0516 | 13   | LS.131.1300 | 49   |
| ED.149.1024        | 31        | FR.073.3818 | 9    | FR.165.3814 | 3    | HK.030.1012 | 13   | LS.151.0315 | 50   |
| ED.150.1024        | 31        | FR.074.3807 | 9    | FR.165.4515 | 3    | HK.030.1512 | 13   | LS.151.0425 | 50   |
| ED.151.0723        | 31        | FR.074.3808 | 9    | FR.172.1918 | 9    | HK.030.2012 | 13   | LS.151.0615 | 50   |
| ED.155.0524        | 33        | FR.075.2210 | 10   | FR.172.2218 | 9    | HK.350.0615 | 13   | LS.151.0625 | 50   |
| ED.160.0124        | 33        | FR.075.3810 | 10   | FR.172.3810 | 9    | HK.390.0616 | 13   | LS.151.0930 | 50   |
| ED.160.0224        | 33        | FR.076.1018 | 10   | FR.172.3818 | 9    | HK.410.0513 | 13   | LS.152.0315 | 50   |
| ED.370.0020        | 32        | FR.076.2018 | 10   | FR.173.1918 | 9    | HK.420.0615 | 13   | LS.152.0615 | 50   |
| ED.371.0519        | 32        | FR.077.7018 | 10   | FR.173.2510 | 9    | HK.431.0718 | 13   | LS.152.0625 | 50   |
| ED.373.0517        | 32        | FR.078.1521 | 10   | FR.176.1018 | 10   | HK.432.0718 | 13   | LS.153.0415 | 50   |
| ED.374.0517        | 32        | FR.079.1519 | 10   | FR.177.7018 | 10   | <b>KN</b>   |      | LS.153.0625 | 50   |
|                    |           | FR.080.1918 | 4    | FR.180.1918 | 4    |             |      | LS.153.0630 | 50   |
|                    |           |             |      |             |      |             |      | LS.160.0215 | 46   |
|                    |           |             |      |             |      |             |      | LS.160.0315 | 46   |



## Index Numerical

| Item        | Page | Item        | Page | Item        | Page | Item        | Page | Item         | Page |
|-------------|------|-------------|------|-------------|------|-------------|------|--------------|------|
| LS.160.0325 | 46   | MA.125.4500 | 34   | SC.062.0612 | 18   | SM.011.0615 | 22   | SP.284.0410  | 21   |
| LS.160.0415 | 46   | MA.125.5000 | 34   | SC.070.0411 | 18   | SM.011.0718 | 22   | SP.286.0411  | 21   |
| LS.160.0225 | 46   | MA.220.3436 | 34   | SC.071.0411 | 18   | SM.021.0411 | 22   | SP.300.0820  | 23   |
| LS.160.0525 | 46   | MA.220.3842 | 34   | SC.072.0411 | 18   | SM.031.0411 | 22   | SP.310.0821  | 23   |
| LS.160.0530 | 46   | MA.220.4550 | 34   | SC.073.0411 | 18   | <b>SP</b>   |      | <b>ST</b>    |      |
| LS.162.0330 | 46   | MA.225.3436 | 34   | SC.080.0411 | 18   | SP.010.3512 | 43   | ST.120.1128  | 44   |
| LS.164.1705 | 47   | MA.225.3842 | 34   | SC.081.0411 | 18   | SP.010.5012 | 43   | ST.225.0615  | 44   |
| LS.164.2005 | 47   | MA.225.4550 | 34   | SC.082.0411 | 18   | SP.010.7512 | 43   | ST.230.0615  | 44   |
| LS.164.0305 | 47   | MA.410.0400 | 34   | SC.083.0411 | 18   | SP.010.9012 | 43   | ST.240.0615  | 44   |
| LS.164.0210 | 47   | MA.955.0014 | 35   | SC.090.0718 | 19   | SP.032.0719 | 17   | ST.325.0615  | 44   |
| LS.164.1715 | 47   | MA.956.0016 | 35   | SC.091.0718 | 19   | SP.090.0410 | 22   | ST.330.0615  | 44   |
| LS.164.0215 | 47   | MA.960.0016 | 35   | SC.095.0410 | 19   | SP.290.0412 | 22   | ST.340.0615  | 44   |
| LS.164.0315 | 47   | MA.960.0116 | 35   | SC.100.0514 | 19   | SP.110.0719 | 16   | <b>TF</b>    |      |
| LS.165.1705 | 47   | MA.965.0014 | 35   | SC.100.0718 | 19   | SP.110.0923 | 16   | TF.010.0412  | 38   |
| LS.165.2005 | 47   | MA.965.0114 | 35   | SC.100.0923 | 19   | SP.110.3512 | 43   | TF.010.0615  | 38   |
| LS.165.0305 | 47   | <b>NS</b>   |      | SC.101.0514 | 19   | SP.110.5012 | 43   | TF.020.0412  | 38   |
| LS.165.0210 | 47   | NS.020.0000 | 43   | SC.101.0718 | 19   | SP.110.7512 | 43   | TF.030.0412  | 38   |
| LS.165.1715 | 47   | NS.025.0023 | 43   | SC.101.0923 | 19   | SP.110.9012 | 43   | TF.032.0412  | 38   |
| LS.165.0215 | 47   | NS.030.0023 | 43   | SC.110.0514 | 19   | SP.111.0719 | 16   | TF.030.0615  | 38   |
| LS.165.0315 | 47   | NS.030.0123 | 43   | SC.110.0718 | 19   | SP.111.0719 | 16   | TF.030.1117  | 41   |
| LS.166.1710 | 47   | NS.030.0223 | 43   | SC.111.0514 | 19   | SP.120.0719 | 16   | TF.030.1223  | 41   |
| LS.166.0210 | 47   | <b>RR</b>   |      | SC.111.0718 | 19   | SP.120.0923 | 16   | TF.030.1728  | 41   |
| LS.166.0215 | 47   | RR.100.0514 | 11   | SC.111.0923 | 19   | SP.121.0719 | 16   | TF.040.0412  | 38   |
| LS.167.0215 | 47   | RR.120.0513 | 11   | SC.120.0514 | 19   | SP.121.0923 | 16   | TF.050.0412  | 38   |
| LS.170.0215 | 47   | RR.130.0513 | 11   | SC.120.0718 | 19   | SP.130.0718 | 16   | TF.050.0514  | 41   |
| LS.170.1215 | 47   | RR.131.0513 | 11   | SC.120.0923 | 19   | SP.131.0718 | 16   | TF.050.0615  | 38   |
| LS.170.1515 | 47   | RR.150.0616 | 11   | SC.121.0514 | 19   | SP.140.0718 | 16   | TF.052.0412  | 38   |
| LS.170.1220 | 47   | RR.155.0616 | 12   | SC.121.0718 | 19   | SP.141.0718 | 16   | TF.052.0514  | 41   |
| LS.170.1520 | 47   | RR.160.0813 | 12   | SC.121.0923 | 19   | SP.150.0718 | 16   | TF.056.0412  | 38   |
| LS.170.2020 | 47   | RR.160.1013 | 12   | SC.131.0512 | 19   | SP.151.0718 | 16   | TF.060.0412  | 39   |
| LS.170.1215 | 47   | RR.160.1213 | 12   | SC.133.0512 | 19   | SP.160.0718 | 16   | TF.060.0615  | 39   |
| LS.170.1530 | 47   | RR.160.1413 | 12   | SC.133.0615 | 19   | SP.170.0718 | 16   | TF.060.0616  | 41   |
| LS.170.0230 | 47   | RR.160.1613 | 12   | SC.140.0309 | 20   | SP.171.0718 | 16   | TF.060.0716  | 41   |
| LS.171.0215 | 47   | RR.160.1813 | 12   | SC.140.0410 | 20   | SP.180.0718 | 16   | TF.060.0821  | 41   |
| LS.171.0230 | 47   | RR.162.1213 | 12   | SC.141.0309 | 20   | SP.181.0718 | 16   | TF.070.0412  | 39   |
| LS.171.1215 | 47   | RR.162.1213 | 12   | SC.141.0410 | 20   | SP.211.0411 | 16   | TF.072.0412  | 39   |
| LS.171.1515 | 47   | RR.162.1413 | 12   | SC.144.0411 | 20   | SP.211.0615 | 16   | TF.077.0412  | 39   |
| LS.171.1220 | 47   | RR.162.1613 | 12   | SC.145.0411 | 20   | SP.211.0718 | 16   | TF.080.0412  | 39   |
| LS.171.1520 | 47   | RR.162.1813 | 12   | SC.148.0411 | 20   | SP.211.0923 | 16   | TF.082.0412  | 39   |
| LS.171.2020 | 47   | RR.164.1014 | 12   | SC.150.0309 | 20   | SP.221.0615 | 16   | TF.084.0412  | 39   |
| LS.171.1530 | 47   | RR.164.1214 | 12   | SC.150.0410 | 20   | SP.221.0718 | 16   | TF.086.0412  | 39   |
| LS.172.0320 | 48   | RR.164.1414 | 12   | SC.151.0309 | 20   | SP.221.0923 | 16   | TF.090.0512  | 40   |
| LS.173.0320 | 48   | RR.164.1614 | 12   | SC.151.0410 | 20   | SP.230.0719 | 17   | TF.100.0411  | 40   |
| LS.174.0320 | 48   | RR.166.0816 | 12   | SC.154.0411 | 20   | SP.240.0617 | 17   | TF.177.0412  | 39   |
| LS.174.0420 | 48   | RR.166.1016 | 12   | SC.155.0411 | 20   | SP.240.0820 | 17   | TF.277.0412  | 39   |
| LS.174.0330 | 48   | RR.166.1216 | 12   | SC.158.0411 | 20   | SP.241.0617 | 17   | <b>TM</b>    |      |
| LS.174.0430 | 48   | RR.166.1416 | 12   | SC.159.0411 | 20   | SP.241.0820 | 17   | TM.005.0410  | 40   |
| LS.175.0320 | 48   | RR.166.1616 | 12   | SC.164.0411 | 20   | SP.250.0617 | 17   | TM.007.0410  | 40   |
| LS.175.0420 | 48   | RR.166.1816 | 12   | SC.170.0616 | 20   | SP.250.0820 | 17   | TM.010.0410  | 40   |
| LS.175.0330 | 48   | RR.170.0821 | 12   | SC.180.0513 | 21   | SP.251.0617 | 17   | TM.011.0410  | 40   |
| LS.175.0430 | 48   | RR.181.0616 | 12   | SC.181.0513 | 21   | SP.251.0820 | 17   | TM.028.0308  | 40   |
| <b>MA</b>   |      | RR.182.0616 | 12   | SC.184.0412 | 22   | SP.260.0617 | 17   | TM.025.0308  | 40   |
| MA.030.0014 | 34   | <b>SC</b>   |      | SC.190.0513 | 21   | SP.260.0820 | 17   | TM.025.0312  | 40   |
| MA.030.0016 | 34   | SC.010.0514 | 17   | SC.191.0513 | 21   | SP.261.0617 | 17   | TM.0411.0001 | 41   |
| MA.035.3800 | 34   | SC.020.0514 | 17   | <b>SH</b>   |      | SP.261.0820 | 17   | TM.0411.0003 | 41   |
| MA.035.4200 | 34   | SC.022.0514 | 17   | SH.103.0514 | 37   | SP.273.0411 | 21   | TM.0411.0004 | 41   |
| MA.035.4500 | 34   | SC.035.0514 | 18   | SH.104.0616 | 37   | SP.275.0411 | 21   | TM.0411.0005 | 41   |
| MA.035.5000 | 34   | SC.045.0514 | 18   | SH.203.0514 | 37   | SP.280.0309 | 21   | TM.0411.0007 | 41   |
| MA.050.4070 | 35   | SC.040.0718 | 18   | <b>SM</b>   |      | SP.281.0309 | 21   | TM.0411.0030 | 41   |
| MA.100.0815 | 35   | SC.041.0718 | 18   | SM.010.0615 | 22   | SP.290.0309 | 21   | TM.0411.0045 | 41   |
| MA.100.0820 | 35   | SC.050.0718 | 18   | SM.010.0718 | 22   | SP.291.0309 | 21   | TM.0411.0050 | 41   |
| MA.125.3800 | 34   | SC.060.0616 | 18   |             |      | SP.282.0410 | 21   |              |      |
| MA.125.4200 | 34   |             |      |             |      | SP.283.0410 | 21   |              |      |



## Index Numerical

| Item         | Page | Item | Page | Item | Page | Item | Page | Item | Page |
|--------------|------|------|------|------|------|------|------|------|------|
| TM.0411.0055 | 41   |      |      |      |      |      |      |      |      |
| TM.0411.0090 | 41   |      |      |      |      |      |      |      |      |
| TP           |      |      |      |      |      |      |      |      |      |
| TP.010.0513  | 40   |      |      |      |      |      |      |      |      |
| TP.020.0513  | 40   |      |      |      |      |      |      |      |      |
| TR           |      |      |      |      |      |      |      |      |      |
| TR.120.1608  | 1    |      |      |      |      |      |      |      |      |
| TR.120.2409  | 1    |      |      |      |      |      |      |      |      |
| TR.120.2410  | 1    |      |      |      |      |      |      |      |      |
| TR.120.2510  | 1    |      |      |      |      |      |      |      |      |
| TR.120.2512  | 1    |      |      |      |      |      |      |      |      |
| TR.120.2515  | 1    |      |      |      |      |      |      |      |      |
| TR.120.3009  | 1    |      |      |      |      |      |      |      |      |
| TR.120.3012  | 1    |      |      |      |      |      |      |      |      |
| TR.120.3015  | 1    |      |      |      |      |      |      |      |      |
| TR.120.3016  | 1    |      |      |      |      |      |      |      |      |
| TR.120.3019  | 1    |      |      |      |      |      |      |      |      |
| TR.220.1608  | 1    |      |      |      |      |      |      |      |      |
| TR.220.2409  | 1    |      |      |      |      |      |      |      |      |
| TR.220.2410  | 1    |      |      |      |      |      |      |      |      |
| TR.220.2510  | 1    |      |      |      |      |      |      |      |      |
| TR.220.2512  | 1    |      |      |      |      |      |      |      |      |
| TR.220.2515  | 1    |      |      |      |      |      |      |      |      |
| TR.220.3009  | 1    |      |      |      |      |      |      |      |      |
| TR.220.3012  | 1    |      |      |      |      |      |      |      |      |
| TR.220.3015  | 1    |      |      |      |      |      |      |      |      |
| TR.220.3016  | 1    |      |      |      |      |      |      |      |      |
| TR.220.3019  | 1    |      |      |      |      |      |      |      |      |

# Feedback



**Instrumed**  
Fine Surgical Instruments

## Plastic Surgery

[illegible]



**Instrumed**  
Fine Surgical Instruments

# PLASTIC SURGERY INSTRUMENTS



65-B, Small Industries Estate,  
P.O. Box 615, Sialkot, Pakistan.

Ph : +92 52 3551769

: +92 52 3559127

Fax: +92 52 3550536

: +92 52 3553799

info@instrumed.com.pk  
sales@instrumed.com.pk  
www.instrumed.com.pk



**Instrumed**  
Fine Surgical Instruments